



By comparing retention behavior in normal-phase high-performance liquid chromatography (NP-HPLC) and supercritical fluid chromatography (SFC), as well as evaluating the selectivity of various stationary phases, we attempted to elucidate the chromatographic properties of SFC. Under NP-HPLC (ethanol/n-hexane) and SFC (ethanol/CO₂) conditions, the retention of hydrophobic, hydrophilic, and ionic compounds was compared. In addition, the hydrophobicity, hydrophilicity, π - π interaction, and dispersion force of 18 different stationary phases were evaluated under SFC conditions.

Column: (HPLC conditions) See figures below (SFC conditions) See figures below
 Column size: 4.6mm.I.D.-250mm (SFC conditions) See figures below
 Mobile phase: See figures below
 Flow rate: 1.0 ml/min
 BPR: --
 Temperature: 40 °C
 Detection: UV254nm
 Equipment: Prominence (Shimadzu) Nexera UC (Shimadzu)






Stationary Phase Selectivity Evaluation

Column name	Stationary phase	(1) Retention	(2) Retention	(3) Hydrophobicity	(4) Hydrophobicity	(5) π - π interaction	(6) Polar group	(7) Dispersion force
COSMOSIL SSL-II	none	0.00	1.06	0.93	1.09	1.74	2.25	1.17
COSMOSIL 5C18-MS-II	Dial Group	0.00	1.50	10.56	0.76	2.88	1.84	0.85
COSMOSIL SPY	Pyridinyl group	0.00	1.10	3.70	0.85	2.52	2.35	0.63
COSMOSIL SHP	Hydroxyphenyl group	0.00	1.31	10.33	1.07	2.28	2.06	0.90
COSMOSIL SHLII	Triazole	0.00	1.76	13.22	0.94	2.44	1.95	0.74
COSMOSIL SFC-MS	Phenyl group	0.00	0.19	2.81	1.36	2.51	1.70	1.21
COSMOSIL 5 π NAP	Naphthyl Group	0.00	0.31	2.79	1.45	3.10	2.10	1.34
COSMOSIL SPYE	Phenyl Group	0.00	0.92	3.54	1.43	3.66	2.39	1.45
COSMOSIL SPBr	Pentabromophenyl group	0.04	0.71	3.34	1.45	5.19	2.34	1.29
COSMOSIL SPPE	Monobromophenyl Group	0.00	0.50	4.70	1.23	2.54	2.21	1.33
COSMOSIL SPFP	Pentafluorophenyl group	0.04	0.15	3.01	1.19	2.13	1.74	1.29
COSMOSIL Scholester	Cholesteryl Group	0.14	0.29	2.38	1.37	3.71	1.64	0.79
COSMOSIL 5C18-MS-II	Octadecyl Group (monomeric)	0.17	0.17	0.89	1.51	2.76	1.05	1.04
COSMOSIL 5C18-AR-II	Octadecyl Group (polymeric)	0.16	0.23	1.09	1.26	2.87	1.05	0.93
2-Strindylene (Not for sale)	Phenyl group	0.00	0.91	6.08	1.02	2.30	2.15	0.92
Picolyamine (Not for sale)	2-Picolyamine	0.00	1.80	11.82	0.89	3.00	2.07	0.85
Aminophenone (Not for sale)	1-Aminophenone	0.00	1.24	5.91	0.82	2.57	2.03	0.76
DEA (Not for sale)	Diethylamine	0.00	2.02	13.85	0.78	2.88	1.78	0.44

Mobile phase: Ethanol/CO₂ = 5/95

In SFC, in addition to the retention properties of normal phase HPLC, secondary interactions, such as π - π interaction and dispersion force, are stronger compared to HPLC, \rightarrow SFC has unique separation potential

- Retention 1: Retention factor (K') of No.1
- Retention 2: k' of No.2
- Hydrophilicity: Separation factor (α) of No.3/No.2
- Hydrophobicity: α of No.4/No.2
- π - π interaction: α of No.5/No.2
- Polar group: α of No.6/No.2
- Dispersion force: α of No.7/No.2



Results and discussion

For hydrophilic and ionic compounds, silica (SL) showed similar retention in NP-HPLC and SFC, while the pyrenylethyl (PYE) and pentabromobenzyl (PBr) phases retained the compounds several times more strongly in SFC compared to NP-HPLC. On the other hand, while hydrophobic compounds were not retained by any of the stationary phases in NP-HPLC, C₁₈, PYE, and PBr retained them slightly in SFC. The hydrophobicity parameters, obtained from evaluating the selectivity of the stationary phases, were, unsurprisingly, negatively correlated with hydrophilicity. On the other hand, dispersion force was positively correlated with hydrophobicity, suggesting that dispersion force is involved in the retention of hydrophobic compounds. From the above, we can surmise that retention in SFC is caused by secondary interactions, such as dispersion force and π - π interaction, working more strongly than in HPLC, in addition to the interactions at work in NP-HPLC.