


Due to the lack of retention on ODS columns, oligosaccharides have typically been analyzed on specialty sugar columns. We have successfully analyzed inulins using COSMOSIL Sugar-D, a specialty column for saccharide analysis, and PBr, a unique reversed-phase column for polar and non-polar analytes.

## Introduction

### About Inulin

Inulin (G-F<sub>n</sub>) is a polysaccharide composed of a glucose unit (G) bonded to multiple fructose units (F). It is an indigestible, water-soluble dietary fiber, abundant in the roots of chicory and tubers of Jerusalem artichokes. Recently, oolong tea with inulin has been released as a functional food in Japan on the basis that inulins feed intestinal bacteria by forming short-chain fatty acids and facilitate the peristaltic movement in the bowels. In the literature as well, there are references to inulin's effect on bowel movements.<sup>(1),(2)</sup>

- (1) Kleessen B. *et al. The American Journal of Clinical Nutrition* **65**(5), 1397-1402(1997).  
(2) Castiglia-Delavaud C. *et al. British Journal of Nutrition* **80**(4), 343-352(1998).



General structure of inulin

## Sample preparation

### Inulin standard

1. Weigh inulin standard.
2. Dissolve in Milli-Q water to a concentration of 10 mg/mL.


Milli-Q® is a trademark of Merck KGaA.

### Oolong tea with inulin

1. Load 400 µL of oolong tea with inulin into a Cosmospin Filter H (pore size 0.45 µm) and centrifuge.
2. Use filtrate as sample.


## Applications

### Analysis using Sugar-D




## Applications

### ● Reversed-phase analysis



### ● Oolong tea with inulin: reversed-phase analysis



## Ordering Information

Application	Product Name	Product grade	Product No.	Size
Columns	COSMOSIL Sugar-D Packed Column	SP (for HPLC)	05397-51	4.6mmI.D.-250mm
	COSMOSIL PBr Packed Column		12395-51	4.6mmI.D.-250mm
	COSMOSIL 5C <sub>18</sub> -PAQ Packed Column		02485-81	4.6mmI.D.-250mm
Filtration	Cosmospin Filter H (0.45 μm)	—	06540-34	100 pkg

For research use only, not intended for diagnostic or drug use.