


# Separating Vitamin D<sub>2</sub> and D<sub>3</sub>, their Metabolites and Epimers



Ken Tseng<sup>1</sup> (ken@nacalaiusa.com), Toshi Ono<sup>1</sup>, Tsunehisa Hirose<sup>2</sup>


<sup>1</sup>Nacalai USA, Inc., San Diego, United States, <sup>2</sup>Nacalai Tesque, Inc., Kyoto, Japan



## Introduction

The accuracy of current vitamin D measurements by immunoassays and LCMS have been questioned due to the overlapping LC peaks with identical m/z values epimers. To solve this problem, we have developed a new HPLC method to achieve baseline separation of vitamin D<sub>2</sub>/D<sub>3</sub>, their 25-OH metabolites and C3-epimers in one single run. A novel core-shell type reversed-phase HPLC column with Cholesterol as the functional group (Cosmocore Cholester) is used in this study. The Cosmocore Cholester column has similar hydrophobicity as regular C18 columns, so it is used in the same exact way. The baseline separation is so complete that it can be used for quantification by UV detector alone at low concentrations. Gradient condition can be employed to further separate vitamin D<sub>2</sub> and vitamin D<sub>3</sub> and their four metabolites/epimers, all in one LCMS (or UV) run.

## Vitamin D Structures




## Vitamin D Metabolite and Epimer Separation – Flow Rate Comparison


Column: COSMOCORE Cholester  
 Column size: 2.1x150mm, 2.6 µm core-shell particles  
 Flow rate: 0.4 mL/min  
 Temperature: 50 °C  
 Detection: UV 265 nm  
 Isocratic: 20A:80B A: 0.1% formic acid in H<sub>2</sub>O B: methanol

### Samples

- 25(OH)D<sub>3</sub>
- epi-25(OH)D<sub>3</sub>
- 25(OH)D<sub>2</sub>
- epi-25(OH)D<sub>2</sub>



## Vitamin D Separation




## Vitamin D, Metabolite, and Epimer Separation


Column: COSMOCORE Cholester  
 Column size: 2.1x150mm, 2.6 µm core-shell particles  
 Flow rate: 0.3 mL/min  
 Temperature: 50 °C  
 Detection: UV 265 nm  
 Data: blank subtraction performed  
 Gradient: A: 0.1% formic acid in H<sub>2</sub>O B: methanol

### Samples

- 25(OH)D<sub>3</sub>
- epi-25(OH)D<sub>3</sub>
- 25(OH)D<sub>2</sub>
- epi-25(OH)D<sub>2</sub>
- vitamin D<sub>2</sub>
- vitamin D<sub>3</sub>



## Vitamin D Metabolite and Epimer LCMS Data



## Conclusion

- Vitamin D<sub>2</sub> and D<sub>3</sub> isocratic separation under 3 minutes using 100% MeOH
- 25(OH) Vitamin D<sub>2</sub> and D<sub>3</sub> metabolites and C-3 epimers were baseline separated under isocratic condition
- All six vitamin D and associated metabolites were separated in a single HPLC gradient run