



## CosmoPAGE Bis-tris Precast Gel

- Bis-tris gel, 15 wells -

12 %, Cat. No. NU001215  
4-12 %, Cat. No. NU041215

### Introduction

CosmoPAGE Bis-Tris Precast Gel is a high-performance and easy to use precast polyacrylamide gel for electrophoresis in Bis-Tris buffer system (MOPS or MES). The optimized gel formula allows CosmoPAGE Bis-Tris Precast Gel to show improved resolution, accurate results, and an extended shelf-life over conventional gel. CosmoPAGE Bis-Tris Precast Gels are available in gradient (4 to 12%) and fixed (12%) concentrations of polyacrylamide in 15-well formats.

### Product Information

Product	Cat. No.	Gel % / # of well	QTY
CosmoPAGE Bis-tris Precast Gel	NU001215	12% / 15 wells	10 gels
CosmoPAGE Bis-tris Precast Gel	NU041215	4 - 12 % / 15 wells	10 gels

### Storage

4 °C

### Shelf Life

12 months

### Features

- Compatible with Bio-Rad Mini-PROTEAN® CORE
- No combs and tape
- Sharp, straight bands
- Numbered well; loading volume up to 28 µl/wells

### Procedures

#### Running Buffer

Dissolve CosmoPAGE MOPS Running Buffer in 1000 ml of deionized water to make the 1X MOPS running buffer. Or prepare 10X stock running buffer with the following recipe:

10X MOPS running buffer:

Tris base 60.6 g, MOPS 104.6 g, SDS 10.0 g, EDTA 3.0 g.  
Deionized water to 1000 ml.

10X MES running buffer preparation:

Tris base 60.6 g, MES 97.6 g, SDS 10.0 g, EDTA 3.0 g.

Deionized water to 1000 ml

\*Always use fresh 1X running buffer

### Sample Preparation

Mix sample with 2X sample buffer. Heat the diluted samples at 90-95°C for 5 min.

Recipe of 2X sample buffer with reducing agents:

62.5 mM Tris-HCl pH 6.8, 2% SDS, 25% (v/v) glycerol, 0.01% bromophenol blue, and 100 mM DTT (or 5% β-mercaptoethanol) as reducing agent.

### Setting CosmoPAGE Precast Gel into Bio-Rad Mini-PROTEAN® Core


Totally the same way to place CosmoPAGE precast gels and Bio-Rad Mini-PROTEAN gels; no adaptors required.

### Running Conditions

	100 V	150 V	200 V
Running time	90 min.	60 min.	35 min.
Current			
Initial per gel	15-20 mA	50-60 mA	100-110 mA
Final per gel	5-10 mA	20-30 mA	30-40 mA
TEMP.	25°C	25-35°C	25-35°C

Running time might vary depending on power supply and the temperature of gel, running buffer and lab.

### Gel Migration chart



Migration patterns of Unstained Protein Markers (#29458-24)

### Gel staining

1. Insert a cassette opener into four corners of cassettes.
2. Carefully **pry but not twist** the opener to separate two plates.
3. Gently remove the gel from the plate.
4. Follow a procedure of most popular staining methods, such as Coomassie dyes (R-250 or G-250), Silver-stain solution, and fluorescent.

### Gel blotting

Follow the general guidelines for your blotting apparatus; compatible to wet, semi-dry, or dry transferred to PVDF or nitrocellulose membrane.

### Technical Support

E-mail: [info@nacalaiusa.com](mailto:info@nacalaiusa.com)

Website: [www.nacalaiusa.com](http://www.nacalaiusa.com)

Rev. 03/2020