

<商標について>  
• COSMOSIL®、コスモシール® は弊社の登録商標です。

● Web site  
<http://www.nacalai.co.jp>

● 価格・納期のご照会  
0120-489-552

● 製品に関する技術的なご照会  
E-mail: info-tech@nacalai.co.jp  
Tel: 075-211-2703 Fax: 075-211-2673

\*希望納入価格に消費税は含まれておりません。  
\*試験・研究以外の用途には使用しないでください。  
\*外観および仕様は改良のため、予告なく変更する場合がございますのでご了承ください。

## ナカライテスク株式会社

本社	604-0855	京都市中京区二条通烏丸西入東玉屋町498	TEL(075)211-2516
本社営業所	604-0855	京都市中京区二条通烏丸西入東玉屋町498	TEL(075)231-5301
京都1営業所	601-8204	京都市南区久世東土川町272	TEL(075)933-6492
京都2営業所	601-8204	京都市南区久世東土川町272	TEL(075)933-6492
京都3営業所	601-8204	京都市南区久世東土川町272	TEL(075)933-6492
滋賀営業所	524-0033	滋賀県守山市浮気町346	TEL(077)581-3610
京阪奈営業所	630-8124	奈良市三条松町5-16	TEL(0742)35-2001
大阪営業所	564-0015	大阪府吹田市幸町23-1	TEL(06)6381-3621
神戸営業所	657-0843	神戸市灘区大石北町7-3	TEL(078)861-0145
福岡営業所	812-0013	福岡市博多区博多駅東3-12-1 アバンダント95ビル (上記営業所窓口)	TEL(092)474-2341
西日本センター	604-0855	京都市中京区二条通烏丸西入東玉屋町498	TEL(075)254-0070

## ■販売取扱店

東京1営業所	103-0021	東京都中央区日本橋本石町4-2-16日本橋本石町トーセイビル	TEL(03)3242-5881
東京2営業所	103-0021	東京都中央区日本橋本石町4-2-16日本橋本石町トーセイビル	TEL(03)3242-5886
神奈川営業所	222-0033	横浜市港北区新横浜1-3-1 新横浜アーバンスクエア	TEL(045)478-2100
埼玉営業所	336-0022	さいたま市南区白鳩1-17-12	TEL(048)866-5851
仙台営業所	981-0915	仙台市青葉区通町2-2-8モンレーブ通町101	TEL(022)727-8070

## HPLC COSMOSIL Application Data HILIC 2012

親水性相互作用クロマトグラフィーは比較的新しい分離機構であり、まだ使用例が多くありません。そこで弊社ではコスモシル HILIC カラムを身近に使っていただくために、分析例をご用意しました。分析したい化合物はもちろんのこと、類似化合物の分離条件の参考としてご利用ください。新しい分析例は随時追加し Web で公開しております。“COSMOSIL アプリケーション” で検索してください。

## CONTENTS

- |                                    |                                       |
|------------------------------------|---------------------------------------|
| 1. 親水性相互作用クロマトグラフィーについて・・・ P2      | 6. COSMOSIL Chromatogram Index・・・ P12 |
| 2. COSMOSIL HILIC・・・ P3            | 7. コスモシル使用者様の声"こんなに役立つ情報"・・・ P33      |
| 3. 移動相の選択について・・・ P6                | 8. Technical Note・・・ P34              |
| 4. COSMOSIL Application Data・・・ P8 | 9. COSMOSIL HILIC掲載文献について・・・ P36      |
| 5. Web Application 検索ページの紹介・・・ P11 | 10. INDEX・・・ P38                      |

### 親水性相互作用クロマトグラフィーについて

親水性相互作用クロマトグラフィー（Hydrophilic Interaction Chromatography=HILIC）は、順相クロマトグラフィーの一種です。逆相クロマトグラフィーとは逆に、カラムに注入されたサンプルは、親水性の低い（疎水性の高い）化合物から溶出します。

分離モード	親水性相互作用クロマトグラフィー	逆相クロマトグラフィー
固定相	親水性基（もしくはシリカゲル）	C18 などの疎水性基
移動相	CH <sub>3</sub> CN などの有機溶媒 / H <sub>2</sub> O	
保持のメカニズム	親水性相互作用	疎水性相互作用
分離対象の化合物	親水性化合物	親水性～疎水性化合物
特長	<ul style="list-style-type: none"> <li>親水性化合物に最適</li> <li>LC-MS に適する</li> </ul>	<ul style="list-style-type: none"> <li>汎用性が高い</li> <li>高分離能</li> <li>分析例（アプリケーション）が豊富</li> </ul>

### 親水性相互作用クロマトグラフィーと逆相クロマトグラフィーの比較


C18 カラムで保持のない親水性化合物を分析する場合、イオンペア試薬を使用することが多いですが、コスモシル HILIC ではイオンペア試薬を使用せずに分離が可能となります。イオンペア試薬の問題点であるカラムの平衡に時間がかかる、カラムの劣化が早い、移動相の調製が手間であるなどが解決できます。また、イオンペア試薬を使わず、揮発性の高い有機溶媒含量が高い移動相を使用することから LC-MS 分析にも適した分離機構です。

### 親水相互作用クロマトグラフィーと逆相クロマトグラフィーとの比較


Column: 4.6mm I.D.-250mm  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: (HILIC) Acetonitrile/10mmol/l Ammonium Acetate = 60/40  
 (C<sub>18</sub>-PAQ) 20mmol/l Phosphate buffer(pH2.5)  
 (Ion pair) 5mmol/l Sodium 7-Octanesulfonate

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1. Glycine (5.0mg/ml)  
 2. Glycylglycine (0.125mg/ml)  
 Inj. Vol 2.0 μl




NACAL TESQUE, INC



# トリアゾール結合型 HPLC カラム COSMOSIL HILIC

- 水素結合 + 陰イオン交換能
- C18 カラムで保持がない親水性化合物の分離に最適
- トリアゾール結合型

固定相構造の一部




トリアゾール


### ● 陰イオン交換能による分離 (他社比較例)

コスモシル HILIC は陰イオン交換能を有するため、酸性化合物に対して保持が大きくなる傾向があります。この分離特性により他社 HILIC カラムとは異なる分離挙動を示します。

陰イオン性に差がある化合物の分離




酸性・塩基性・中性化合物の分離



## ● 親水性相互作用による分離


コスモシル HILIC は固定相のトリアゾールの効果による親水性相互作用（水素結合による）と陰イオン交換能の2つの相互作用によって化合物を分離します。2つの相互作用は移動相によって使い分けられます。

### 陰イオン性に差がない化合物の分離



## ● LC-MS を用いたメラミンの分析例

イオンペア試薬を含まない移動相条件で使用できますので、LC-MS にも適用可能です。




上記 LC-MS/MS データはお客様よりご提供いただきました。

## ● 応用例

C<sub>18</sub> カラムと組み合わせるとお使いいただくことにより分離の幅が広がります。

### C<sub>18</sub> との併用利用



## ● 価格表

### COSMOSIL HILIC Packed Column

サイズ 内径×長さ(mm)	製品番号	価格
1.0 × 150	07869-11	70,000
1.0 × 250	07870-71	82,000
2.0 × 30	08568-21	45,000
2.0 × 50	07052-91	45,000
2.0 × 100	08569-11	50,000
2.0 × 150	07054-71	60,000
2.0 × 250	07489-91	72,000
3.0 × 150	07871-61	50,000
3.0 × 250	07872-51	63,000

### COSMOSIL HILIC Guard Column

サイズ 内径×長さ(mm)	製品番号	価格
4.6 × 150	07056-51	50,000
4.6 × 150 3口セット	09385-23	150,000
4.6 × 250	07057-41	63,000
10 × 250	07059-21	150,000
20 × 250	07060-81	350,000
28 × 250	07875-21	ご照会
2.0 × 10	07569-41	45,000
4.6 × 10	07055-61	25,000
10 × 20	07058-31	53,000
20 × 20	07854-91	115,000
20 × 50	07873-41	145,000
28 × 50	07874-31	ご照会

各種カラムサイズを揃えていますので、詳しくはご照会ください。

## 移動相の選択について

コスモシル HILIC カラムは、親水性相互作用（おもに水素結合による）と陰イオン交換能によって、保持と分離を行います。下記をご覧ください、適切な移動相条件を選択ください。

### (1) 有機溶媒の種類と濃度の効果

- 一般には、アセトニトリル / 水系の移動相を使用します。
- アセトニトリル濃度を高くすると保持は大きくなり、低くすると保持は小さくなります。(図1)
- アセトニトリル濃度は0 ~ 95% (通常50 ~ 95%) の範囲でご使用ください。
- メタノール / 水系は保持が小さくなります。(図2)
- 有機溶媒および水は必ず高速液体クロマトグラフ用をご使用ください。




図1. アセトニトリル濃度と保持との関係  
Column: COSMOSIL HILIC  
Mobile phase: Acetonitrile / 10mmol/l CH<sub>3</sub>COONH<sub>4</sub>




図2. アセトニトリルとメタノールの保持力の違い  
Column: COSMOSIL HILIC  
Mobile phase: Organic solvent / 10mmol/l CH<sub>3</sub>COONH<sub>4</sub> = 90/10

### (2) pH の効果

- 移動相の pH は 2 ~ 7.5 の範囲でご使用ください。
- 中性付近の pH で使用すると保持は大きくなります。(図3)




図3. pH と保持との関係  
Column: COSMOSIL HILIC  
Mobile phase: Acetonitrile / 10mmol/l bu er = 90/10

### (3) 緩衝液の種類と濃度の効果

- 分離性の試料を分析する場合には、移動相に塩や緩衝液を添加する必要があります。
- HILIC では高濃度のアセトニトリルを移動相としますが、アセトニトリルは塩に対する溶解性が低いので注意が必要です。逆相クロマトグラフィーで汎用されるりん酸塩は溶解性が低いために HILIC で使用される場合は注意してください。もし、使用される場合には、アセトニトリル濃度を 70% 以下にしてください。また、ご使用前に移動相に塩が析出していないか、必ず確認してください。
- 緩衝液の種類としては、高濃度のアセトニトリル中でも溶解性の高い、酢酸アンモニウムやギ酸アンモニウムを推奨します。
- 緩衝液の濃度は、5 ~ 100mmol/l の範囲でご使用ください。また、アセトニトリルと混合後、塩の析出がないか確認してください。
- 塩濃度を高くすると、イオン交換能が抑制され、保持が小さくなります。(図4)
- 移動相に緩衝液を使用する場合には、使用前に必ず 0.45 μ m 以下のフィルターでろ過してください。




図4. 塩濃度と保持との関係  
Column: COSMOSIL HILIC  
Mobile phase: Acetonitrile / CH<sub>3</sub>COONH<sub>4</sub> = 50/50

### (4) 移動相選択フロー

ここでは、化合物の特性から、第一選択の移動相を提示致します。保持時間の増減は主にアセトニトリル濃度で調整してください。


・中性化合物	Acetonitrile / Water = 90/10
・塩基性化合物	Acetonitrile / 10mmol/l CH <sub>3</sub> COONH <sub>4</sub> = 90/10
・両性化合物	Acetonitrile / 10mmol/l CH <sub>3</sub> COONH <sub>4</sub> = 70/30
・酸性化合物	Acetonitrile / 10mmol/l CH <sub>3</sub> COONH <sub>4</sub> = 50/50
	溶ししない場合
	Acetonitrile / 10mmol/l Phosphate Bu er(pH7.0) = 50/50

### (5) 分離機能の使い分け

コスモシル HILIC カラムは、2つの分離機能（親水性相互作用と陰イオン交換）によって、保持と分離を行います。移動相を選択することにより、分離機能を使い分けすることができます。陰イオン交換能は塩濃度によって調整できます。塩濃度を上げることで陰イオン交換能を抑制することが可能です。一方、親水性相互作用は、アセトニトリル濃度によって調整できます。アセトニトリル濃度を上げることで親水性相互作用を強めることができます。下図では、酸性化合物であるアスコルビン酸とイソアスコルビン酸の分離について示しています。(4)の項で推奨している移動相 (Acetonitrile / 10mmol/l CH<sub>3</sub>COONH<sub>4</sub> = 50/50) では、陰イオン交換能が強く働き保持が生じていますが、2つの化合物間の陰イオン性の差が少ないため、分離が不十分です。そのため、塩濃度を高めて陰イオン交換能を抑え、さらに親水性相互作用を強めるためにアセトニトリル濃度を上げました。その結果、親水性相互作用によって十分な分離が達成できました。

#### 親水性相互作用による分離

50%Acetonitrile / 10mM CH<sub>3</sub>COONH<sub>4</sub>      50%Acetonitrile / 100mM CH<sub>3</sub>COONH<sub>4</sub>      80%Acetonitrile / 100mM CH<sub>3</sub>COONH<sub>4</sub>  
塩濃度を高めて陰イオン交換能を抑制



### (6) テーリング時の対処法

ピークがテーリングする場合には、下記の移動相を試してください。ピーク形状が改善する場合があります。

- 移動相に 5mmol/l の EDTA を添加する。
  - くえん酸緩衝液にする。(例 10mmol/l Citrate bu er(pH7.0))
- 下記では 5mmol/l EDTA を添加した例を示します。




図5. ピーク形状の改善  
Column: COSMOSIL HILIC(4.6mm I.D.-250mm)  
Flow rate: 1.0ml/min  
Temperature: 30  
Detection: UV254nm  
Sample: Tryptophan(1ng)


### (7) カラムの洗浄

ベースラインが乱れる場合には、50% アセトニトリルを約 30 分間送液して洗浄してください。

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium Acetate = 70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV225nm

Sample:  
1; Melamine (0.1 µg)  
2; Ammeline (0.075 µg)  
3; Cyanuric Acid (0.75 µg)  
4; Ammelide (0.05 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220nm

Sample:  
3; Angiotensin II, Des-Asp1-[Ile<sup>6</sup>] (0.5 µg)  
4; Angiotensin II, [Sar<sup>1</sup>,Ile<sup>2</sup>] (0.5 µg)  
5; Angiotensin II, [Asn<sup>1</sup>,Val<sup>2</sup>] (0.5 µg)  
6; Angiotensin II, [Val<sup>2</sup>] (0.5 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 50mmol/l Ammonium Acetate = 80/20  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD

Sample:  
1; meso-Erythritol (3.0 µg)  
2; Tris(hydroxymethyl)aminomethane (3.0 µg)  
3; Glyceric Acid (6.0 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample:  
1; Trimethylene Glycol (20 µg)  
2; Ethylene Glycol (20 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample:  
1; Ascorbic Acid (1.5 µg)  
2; Malic Acid (3.0 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220nm

Sample:  
1; Nicotinamide (0.125 µg)  
2; Pyridoxine(Vitamin B<sub>6</sub>) (0.25 µg)  
3; Riboflavin (Vitamin B<sub>2</sub>) (0.25 µg)  
4; Nicotinic Acid (0.125 µg)  
5; D-Pantothenic Acid (3.125 µg)  
6; L(+)-Ascorbic Acid (0.875 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7.0) = 60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample:  
1; D-Fructose-6-phosphate (50 µg)  
2; D-Glucose-6-phosphate (50 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7.0) = 60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample:  
1; Glucose (10 µg)  
2; α-D-Glucose-1-phosphate (50 µg)  
3; D-Glucose-6-phosphate (50 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample:  
1; Sorbic Acid (0.05 µg)  
2; Isoascorbic Acid (0.3 µg)  
3; Ascorbic Acid (0.3 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 2.0mmI.D.-150mm  
Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
Flow rate: 0.2 ml/min  
Temperature: 30°C  
Detection: UV220nm

Sample:  
1; Nicotinamide (0.125 µg)  
2; Pyridoxine(Vitamin B<sub>6</sub>) (0.25 µg)  
3; Riboflavin (Vitamin B<sub>2</sub>) (0.25 µg)  
4; Nicotinic Acid (0.125 µg)  
5; D-Pantothenic Acid (3.125 µg)  
6; L(+)-Ascorbic Acid (0.875 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample:  
1; Trimethylene Glycol (20 µg)  
2; Glycerol (20 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium Acetate = 60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample:  
1; Glycine (10 µg)  
2; Glycylglycine (0.25 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7.0) = 70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample:  
1; L-Citrulline (7.5 µg)  
2; Malic Acid (3.8 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile / H<sub>2</sub>O = 95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI

Sample:  
1; Diethylene Glycol (20 µg)  
2; Glycerol (20 µg)




NACALAI TESQUE, INC

### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium Acetate = 85/15  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD

Sample:  
1; Leucine (3.0 µg)  
2; Isoleucine (3.0 µg)  
3; Valine (3.0 µg)




NACALAI TESQUE, INC













### COSMOSIL Application Data

Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210nm

Sample:  
1; Oxamic Acid (0.2 µg)  
2; Oxalic Acid (1.0 µg)



NACALAI TESQUE, INC



**COSMOSIL Application** (<http://www.nacalai.co.jp/cosmosil/data/csmsrchttop.cfm>)

■分析例：業界最大級の7000件以上  
 分析したい化合物の分析は、どのカラム・どんな分析条件で行えばいいの？と迷った際、豊富な分析例の中からの確かなカラム・分析条件を見つけることができます。“COSMOSIL アプリケーション”で検索してください。

■簡単検索  
 サンプル名は、もちろん、分析対象のカテゴリー、カラム名、CAS 番号、粒子径で検索が可能です。化合物のいろいろな情報からお探しの分析例を簡単に見つけることができます。

検索画面

検索結果画面

各データ


●複数のサンプルの同時分析例や異なるカラムでの比較例などを掲載

●化合物の個々の分析例・構造を掲載

今回 COSMOSIL Application では、従来の Application data (A#-\*\*\*\*) と COSMOSIL Chromatogram INDEX (C<-\*\*\*\*) を統合いたしました。


### COSMOSIL Chromatogram Index

Sample: Acetamide  
CAS No.: [60-35-5]  
Molecular formula: C<sub>2</sub>H<sub>5</sub>NO  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O=95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 4.75min  
Capacity factor: 0.57




### COSMOSIL Chromatogram Index

Sample: Acetazolamide  
CAS No.: [59-66-5]  
Molecular formula: C<sub>4</sub>H<sub>7</sub>N<sub>5</sub>O<sub>2</sub>S<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.2mg/ml  
Injection volume: 5.90µl  
Retention time: 5.99min  
Capacity factor: 1.05




### COSMOSIL Chromatogram Index

Sample: p-Aminobenzoic Acid  
CAS No.: [150-13-0]  
Molecular formula: C<sub>7</sub>H<sub>7</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.4mg/ml  
Injection volume: 1.0µl  
Retention time: 16.97min  
Capacity factor: 4.91




### COSMOSIL Chromatogram Index

Sample: 4-Amino-n-butyric Acid  
CAS No.: [56-12-2]  
Molecular formula: C<sub>4</sub>H<sub>7</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 7.67min  
Capacity factor: 1.92




### COSMOSIL Chromatogram Index

Sample: Acetizic Acid  
CAS No.: [85-36-9]  
Molecular formula: C<sub>8</sub>H<sub>7</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.8mg/ml  
Injection volume: 1.0µl  
Retention time: 16.39min  
Capacity factor: 4.76




### COSMOSIL Chromatogram Index

Sample: Acrylic Acid  
CAS No.: [79-10-7]  
Molecular formula: C<sub>3</sub>H<sub>4</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 15.05min  
Capacity factor: 4.28




### COSMOSIL Chromatogram Index

Sample: 6-Amino-hexanoic Acid  
CAS No.: [60-32-2]  
Molecular formula: C<sub>6</sub>H<sub>11</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 7.98min  
Capacity factor: 2.05




### COSMOSIL Chromatogram Index

Sample: 5-Amino-levulinic Acid  
CAS No.: [5451-09-2]  
Molecular formula: C<sub>5</sub>H<sub>7</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0µl  
Retention time: 12.80min  
Capacity factor: 3.87




### COSMOSIL Chromatogram Index

Sample: L-α-Alanine  
CAS No.: [56-41-7]  
Molecular formula: C<sub>3</sub>H<sub>7</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 2.0µl  
Retention time: 6.67min  
Capacity factor: 1.53




### COSMOSIL Chromatogram Index

Sample: β-Alanine  
CAS No.: [107-95-9]  
Molecular formula: C<sub>3</sub>H<sub>5</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 7.38min  
Capacity factor: 1.81




### COSMOSIL Chromatogram Index

Sample: 2-Aminopyridine  
CAS No.: [504-29-0]  
Molecular formula: C<sub>5</sub>H<sub>5</sub>N<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 4.25min  
Capacity factor: 0.39




### COSMOSIL Chromatogram Index

Sample: 3-Aminopyridine  
CAS No.: [462-08-8]  
Molecular formula: C<sub>5</sub>H<sub>5</sub>N<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 1.0µl  
Retention time: 4.05min  
Capacity factor: 0.51




### COSMOSIL Chromatogram Index

Sample: Allantoin  
CAS No.: [97-59-6]  
Molecular formula: C<sub>4</sub>H<sub>6</sub>N<sub>2</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 8.75min  
Capacity factor: 2.02




### COSMOSIL Chromatogram Index

Sample: p-Aminobenzamide  
CAS No.: [3858-83-1]  
Molecular formula: C<sub>7</sub>H<sub>8</sub>N<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Citrate buffer(pH7.0)=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 6.07min  
Capacity factor: 1.31




### COSMOSIL Chromatogram Index

















Sample: 5-Amino-1H-tetrazole  
CAS No.: [4418-61-5]  
Molecular formula: CH<sub>2</sub>N<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 15.76min  
Capacity factor: 4.49



### COSMOSIL Chromatogram Index

Sample: 3-Amino-1H-1,2,4-triazole  
CAS No.: [61-82-5]  
Molecular formula: C<sub>2</sub>H<sub>3</sub>N<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.2mg/ml  
Injection volume: 1.0µl  
Retention time: 7.01min  
Capacity factor: 1.42









### COSMOSIL Chromatogram Index

Sample: Benzamidine  
CAS No.: [618-39-3]  
Molecular formula:  $C_8H_{10}N_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.5mg/ml  
Injection volume: 0.5µl  
Retention time: 7.16min  
Capacity factor: 1.46




### COSMOSIL Chromatogram Index

Sample: Benzenesulfonic Acid  
CAS No.: [98-11-3]  
Molecular formula:  $C_6H_5O_3S$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0µl  
Retention time: 11.54min  
Capacity factor: 3.05




### COSMOSIL Chromatogram Index

Sample: Ceftriaxone  
CAS No.: [73384-59-5]  
Molecular formula:  $C_{18}H_{16}N_4O_5S_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 11.36min  
Capacity factor: 3.05




### COSMOSIL Chromatogram Index

Sample: Chloroacetic Acid  
CAS No.: [79-11-8]  
Molecular formula:  $C_2H_3ClO_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 14.69min  
Capacity factor: 4.15




### COSMOSIL Chromatogram Index

Sample: Benzoic Acid  
CAS No.: [65-85-0]  
Molecular formula:  $C_7H_6O_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 5.0mg/ml  
Injection volume: 0.5µl  
Retention time: 15.19min  
Capacity factor: 4.29




### COSMOSIL Chromatogram Index

Sample: Benzylamine  
CAS No.: [100-46-9]  
Molecular formula:  $C_8H_{11}N$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/50mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 5.58min  
Capacity factor: 0.95




### COSMOSIL Chromatogram Index

Sample: Citrazinic Acid  
CAS No.: [99-11-6]  
Molecular formula:  $C_8H_6NO_4$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 14.16min  
Capacity factor: 3.98




### COSMOSIL Chromatogram Index

Sample: Creatine  
CAS No.: [57-00-1]  
Molecular formula:  $C_4H_9N_3O_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 6.35min  
Capacity factor: 1.40




### COSMOSIL Chromatogram Index

Sample: Bromoacetic Acid  
CAS No.: [79-08-3]  
Molecular formula:  $C_2H_3BrO_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 15.04min  
Capacity factor: 4.31




### COSMOSIL Chromatogram Index

Sample: Cacotheline  
CAS No.: [561-20-6]  
Molecular formula:  $C_{12}H_{17}N_3O_7$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 12.19min  
Capacity factor: 3.23




### COSMOSIL Chromatogram Index

Sample: Creatinine  
CAS No.: [60-27-5]  
Molecular formula:  $C_4H_7N_3O$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 8.93min  
Capacity factor: 2.08




### COSMOSIL Chromatogram Index

Sample: Cyanoacetic Acid  
CAS No.: [372-09-8]  
Molecular formula:  $C_3H_3NO_2$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6, Atten=8  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 11.78min  
Capacity factor: 3.56




### COSMOSIL Chromatogram Index

Sample: Camostat  
CAS No.: [59721-28-7]  
Molecular formula:  $C_{22}H_{26}N_2O_5$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV265 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 7.16min  
Capacity factor: 1.47




### COSMOSIL Chromatogram Index

Sample: L-Carnitine  
CAS No.: [541-15-1]  
Molecular formula:  $C_7H_{15}NO_3$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6, Atten=8  
Sample conc.: 2.0mg/ml  
Injection volume: 1.5µl  
Retention time: 6.96min  
Capacity factor: 1.78




### COSMOSIL Chromatogram Index

Sample: Cyanuric Acid  
CAS No.: [108-80-5]  
Molecular formula:  $C_3H_3N_3O_3$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 7.61min  
Capacity factor: 1.68




### COSMOSIL Chromatogram Index

Sample: L-Cysteine  
CAS No.: [52-90-4]  
Molecular formula:  $C_3H_7NO_2S$   
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Phosphate buffer(pH7)=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 5.0mg/ml  
Injection volume: 2.0µl  
Retention time: 7.05min  
Capacity factor: 1.69



### COSMOSIL Chromatogram Index


Sample: L-(+)-Cystine  
CAS No.: [56-89-3]  
Molecular formula: C<sub>4</sub>H<sub>8</sub>N<sub>2</sub>O<sub>4</sub>S<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)-60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 0.5µl  
Retention time: 6.42min  
Capacity factor: 1.38



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Cytidine  
CAS No.: [65-46-3]  
Molecular formula: C<sub>8</sub>H<sub>11</sub>N<sub>3</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-80/20  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 7.22min  
Capacity factor: 1.58



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Datrizoic Acid  
CAS No.: [117-96-4]  
Molecular formula: C<sub>7</sub>H<sub>7</sub>N<sub>3</sub>O<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.8mg/ml  
Injection volume: 1.0µl  
Retention time: 14.98min  
Capacity factor: 4.26



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Dipicolinic acid  
CAS No.: [499-83-2]  
Molecular formula: C<sub>5</sub>H<sub>3</sub>N<sub>2</sub>O<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)-50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 6.37min  
Capacity factor: 1.23



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Cytosine  
CAS No.: [71-30-7]  
Molecular formula: C<sub>4</sub>H<sub>5</sub>N<sub>3</sub>O  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 11.22min  
Capacity factor: 2.87



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 3,4-Diaminobenzoic Acid  
CAS No.: [619-05-6]  
Molecular formula: C<sub>7</sub>H<sub>7</sub>N<sub>3</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/20mmol/l Ammonium acetate-50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.10mg/ml  
Injection volume: 4.0µl  
Retention time: 16.13min  
Capacity factor: 4.62



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Dithiouracil  
CAS No.: [2001-93-6]  
Molecular formula: C<sub>4</sub>H<sub>4</sub>N<sub>2</sub>S<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.2mg/ml  
Injection volume: 1.5µl  
Retention time: 15.60min  
Capacity factor: 4.15



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: L-DOPA  
CAS No.: [59-92-7]  
Molecular formula: C<sub>9</sub>H<sub>9</sub>NO<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)-70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 3.0mg/ml  
Injection volume: 3.0µl  
Retention time: 7.12min  
Capacity factor: 1.72



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 3,5-Diaminobenzoic Acid  
CAS No.: [535-87-5]  
Molecular formula: C<sub>7</sub>H<sub>7</sub>N<sub>3</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/20mmol/l Ammonium acetate-50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 4.0µl  
Retention time: 16.54min  
Capacity factor: 4.76



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 2,4-Diaminophenol  
CAS No.: [95-86-3]  
Molecular formula: C<sub>6</sub>H<sub>7</sub>N<sub>2</sub>O  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/20mmol/l Ammonium acetate-90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 4.40min  
Capacity factor: 0.51



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Dopamine  
CAS No.: [51-61-6]  
Molecular formula: C<sub>8</sub>H<sub>9</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)-70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 4.0µl  
Retention time: 7.73min  
Capacity factor: 1.96



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: L-(+)-Ergothioneine  
CAS No.: [497-30-3]  
Molecular formula: C<sub>8</sub>H<sub>11</sub>N<sub>3</sub>O<sub>4</sub>S  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-80/20  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 10.29min  
Capacity factor: 2.79



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: DL-2,6-Diaminopimelic Acid  
CAS No.: [583-93-7]  
Molecular formula: C<sub>8</sub>H<sub>11</sub>N<sub>3</sub>O<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)-60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.5µl  
Retention time: 6.93min  
Capacity factor: 1.56



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: DL-2,3-Diaminopropionic Acid  
CAS No.: [54897-59-5]  
Molecular formula: C<sub>3</sub>H<sub>6</sub>N<sub>2</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6,Atten=8  
Sample conc.: 5.0mg/ml  
Injection volume: 2.0µl  
Retention time: 14.38min  
Capacity factor: 4.52



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: meso-Erythritol [meso-Erythrite]  
CAS No.: [149-32-6]  
Molecular formula: C<sub>4</sub>H<sub>8</sub>O<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/H<sub>2</sub>O-90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6,Atten=8  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.78min  
Capacity factor: 1.18



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Famotidin  
CAS No.: [76824-35-6]  
Molecular formula: C<sub>14</sub>H<sub>14</sub>N<sub>4</sub>O<sub>2</sub>S  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate-90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.25mg/ml  
Injection volume: 2.0µl  
Retention time: 5.99min  
Capacity factor: 1.06



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: Folic Acid  
CAS No.: [59-30-3]  
Molecular formula:  $C_{19}H_{19}N_7O_6$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.25mg/ml  
Injection volume: 2.0µl  
Retention time: 14.09min  
Capacity factor: 3.95




### COSMOSIL Chromatogram Index

Sample: Folinic Acid  
CAS No.: [58-05-9]  
Molecular formula:  $C_{20}H_{21}N_7O_7$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.25mg/ml  
Injection volume: 2.0µl  
Retention time: 13.36min  
Capacity factor: 3.68




### COSMOSIL Chromatogram Index

Sample:  $\alpha$ -D-Glucose-1-phosphate  
CAS No.: [59-56-3]  
Molecular formula:  $C_6H_{12}O_6P$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI  
Attenuation:  $4 \times 10^3$  RIU/FS  
Sample conc.: 10.0mg/ml  
Injection volume: 5.0µl  
Retention time: 12.26min  
Capacity factor: 3.68




### COSMOSIL Chromatogram Index

Sample: D-Gluconic Acid  
CAS No.: [6556-12-3]  
Molecular formula:  $C_6H_{12}O_7$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.45min  
Capacity factor: 0.92




### COSMOSIL Chromatogram Index

Sample: Formamide  
CAS No.: [75-12-7]  
Molecular formula:  $CH_3NO$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O=95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 4.58min  
Capacity factor: 0.52




### COSMOSIL Chromatogram Index

Sample: D-Fructose-6-phosphate  
CAS No.: [643-13-0]  
Molecular formula:  $C_6H_{12}O_6P$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI  
Attenuation:  $4 \times 10^3$  RIU/FS  
Sample conc.: 10.0mg/ml  
Injection volume: 5.0µl  
Retention time: 12.16min  
Capacity factor: 3.64




### COSMOSIL Chromatogram Index

Sample: L-Glutamic Acid  
CAS No.: [56-86-0]  
Molecular formula:  $C_5H_9NO_4$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 19.38min  
Capacity factor: 5.87




### COSMOSIL Chromatogram Index

Sample: L-Glutamine  
CAS No.: [56-85-9]  
Molecular formula:  $C_5H_{11}N_2O_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 7.50min  
Capacity factor: 1.85




### COSMOSIL Chromatogram Index

Sample: Fuchsine, Acid  
CAS No.: [3244-88-0]  
Molecular formula:  $C_{20}H_{11}N_3Na_2O_5S_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 1.5µl  
Retention time: 13.82min  
Capacity factor: 3.85




### COSMOSIL Chromatogram Index

Sample: Fumaric Acid  
CAS No.: [110-17-8]  
Molecular formula:  $C_4H_4O_4$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.2mg/ml  
Injection volume: 0.5µl  
Retention time: 10.63min  
Capacity factor: 2.75




### COSMOSIL Chromatogram Index

Sample: Glutaric Acid  
CAS No.: [110-94-1]  
Molecular formula:  $C_5H_8O_4$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 10.45min  
Capacity factor: 2.68




### COSMOSIL Chromatogram Index

Sample: Glutathione(Reduced Form)  
CAS No.: [70-18-8]  
Molecular formula:  $C_9H_{16}N_2O_6S$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 2.0µl  
Retention time: 5.43min  
Capacity factor: 0.89




### COSMOSIL Chromatogram Index

Sample: Gluconic Acid  
CAS No.: [526-95-4]  
Molecular formula:  $C_6H_{12}O_7$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.15min  
Capacity factor: 0.81




### COSMOSIL Chromatogram Index

Sample: D-Glucose-6-phosphate  
CAS No.: [56-73-5]  
Molecular formula:  $C_6H_{12}O_6P$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI  
Attenuation:  $4 \times 10^3$  RIU/FS  
Sample conc.: 10.0mg/ml  
Injection volume: 5.0µl  
Retention time: 12.95min  
Capacity factor: 3.94




### COSMOSIL Chromatogram Index

Sample: DL-Glyceric Acid  
CAS No.: [600-19-1]  
Molecular formula:  $C_3H_6O_4$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 6.0mg/ml  
Injection volume: 5.0µl  
Retention time: 15.68min  
Capacity factor: 4.50




### COSMOSIL Chromatogram Index

Sample: Glycinamide  
CAS No.: [588-14-1]  
Molecular formula:  $C_2H_5N_2O$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELS  
Attenuation: Gain=6, Atten=8  
Sample conc.: 1.0mg/ml  
Injection volume: 3.0µl  
Retention time: 13.35min  
Capacity factor: 3.64



### COSMOSIL Chromatogram Index


Sample: Glycine  
CAS No.: [56-40-6]  
Molecular formula: C<sub>2</sub>H<sub>5</sub>N<sub>2</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 2.0µl  
Retention time: 7.29min  
Capacity factor: 1.77



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Glycolic Acid  
CAS No.: [79-14-1]  
Molecular formula: C<sub>2</sub>H<sub>3</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 15.28min  
Capacity factor: 4.39



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Hydanotic Acid  
CAS No.: [462-60-2]  
Molecular formula: C<sub>3</sub>H<sub>4</sub>N<sub>2</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0µl  
Retention time: 16.33min  
Capacity factor: 4.72



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Hydanoin  
CAS No.: [461-72-3]  
Molecular formula: C<sub>3</sub>H<sub>4</sub>N<sub>2</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 5.01min  
Capacity factor: 0.66



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Glycylglycine  
CAS No.: [556-50-3]  
Molecular formula: C<sub>4</sub>H<sub>8</sub>N<sub>2</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 6.40min  
Capacity factor: 1.27



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Guanidoacetic Acid  
CAS No.: [352-97-6]  
Molecular formula: C<sub>3</sub>H<sub>6</sub>N<sub>2</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 6.61min  
Capacity factor: 1.51



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Hydroxylamine-*O*-sulfonic Acid  
CAS No.: [2950-43-8]  
Molecular formula: H<sub>3</sub>N<sub>2</sub>O<sub>3</sub>S  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELS  
Attenuation: Gain=6,Atten=8  
Sample conc.: 2.0mg/ml  
Injection volume: 3.0µl  
Retention time: 15.60min  
Capacity factor: 5.24



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *cis*-4-Hydroxy-*D*-proline  
CAS No.: [2584-71-6]  
Molecular formula: C<sub>5</sub>H<sub>8</sub>N<sub>2</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 6.96min  
Capacity factor: 1.65



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 1,2,6-Hexanetriol  
CAS No.: [106-69-4]  
Molecular formula: C<sub>6</sub>H<sub>14</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELS  
Attenuation: Gain=6,Atten=8  
Sample conc.: 1.0mg/ml  
Injection volume: 2.0µl  
Retention time: 5.19min  
Capacity factor: 0.80



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *L*-Histidine  
CAS No.: [71-00-1]  
Molecular formula: C<sub>6</sub>H<sub>9</sub>N<sub>3</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.50mg/ml  
Injection volume: 1.0µl  
Retention time: 8.38min  
Capacity factor: 2.19



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *L*-Hydroxyproline  
CAS No.: [51-35-4]  
Molecular formula: C<sub>5</sub>H<sub>9</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 6.49min  
Capacity factor: 1.47



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *N*-Hydroxysuccinimide  
CAS No.: [6066-82-6]  
Molecular formula: C<sub>4</sub>H<sub>5</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 1.5µl  
Retention time: 11.29min  
Capacity factor: 3.22



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *L*-Homocystine  
CAS No.: [626-72-2]  
Molecular formula: C<sub>8</sub>H<sub>12</sub>N<sub>2</sub>O<sub>6</sub>S<sub>2</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 2.0mg/ml  
Injection volume: 1.0µl  
Retention time: 13.41min  
Capacity factor: 4.10



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *L*-Homoserine  
CAS No.: [672-15-1]  
Molecular formula: C<sub>4</sub>H<sub>9</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 7.03min  
Capacity factor: 1.67



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Indigo carmine  
CAS No.: [860-22-0]  
Molecular formula: C<sub>16</sub>H<sub>8</sub>N<sub>2</sub>Na<sub>2</sub>O<sub>6</sub>S<sub>2</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.2mg/ml  
Injection volume: 1.0µl  
Retention time: 7.82min  
Capacity factor: 1.79



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *D*-Isoscorbic Acid  
CAS No.: [89-65-6]  
Molecular formula: C<sub>6</sub>H<sub>8</sub>O<sub>6</sub>  
Column: HILIC  
Column size: 4.6mmI.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 245nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.2mg/ml  
Injection volume: 3.0µl  
Retention time: 17.26min  
Capacity factor: 5.11



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: Isonichomeronic Acid  
CAS No.: [100-26-5]  
Molecular formula:  $C_6H_7NO_4$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 5.91min  
Capacity factor: 1.07




### COSMOSIL Chromatogram Index

Sample: L-Isoleucine  
CAS No.: [73-32-5]  
Molecular formula:  $C_6H_{13}NO_2$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.12min  
Capacity factor: 0.95




### COSMOSIL Chromatogram Index

Sample: L-Lysine  
CAS No.: [56-87-1]  
Molecular formula:  $C_6H_{13}N_2O_2$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/50mmol/l Phosphate buffer(pH7.0)=60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 9.26min  
Capacity factor: 2.55




### COSMOSIL Chromatogram Index

Sample: Melamine  
CAS No.: [108-78-1]  
Molecular formula:  $C_3H_6N_6$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV240 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 10.79min  
Capacity factor: 2.79




### COSMOSIL Chromatogram Index

Sample: Isonicotinic Acid  
CAS No.: [55-22-1]  
Molecular formula:  $C_6H_5NO_2$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 16.45min  
Capacity factor: 4.78




### COSMOSIL Chromatogram Index

Sample: Isonicotinohydrazide  
CAS No.: [54-85-3]  
Molecular formula:  $C_6H_7N_3O$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV265 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.50mg/ml  
Injection volume: 0.5µl  
Retention time: 5.37min  
Capacity factor: 0.85




### COSMOSIL Chromatogram Index

Sample: Maleic Acid  
CAS No.: [110-16-7]  
Molecular formula:  $C_4H_2O_4$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 7.45min  
Capacity factor: 1.62




### COSMOSIL Chromatogram Index

Sample: L-(+)-Malic Acid  
CAS No.: [97-67-6]  
Molecular formula:  $C_4H_6O_5$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 10.55min  
Capacity factor: 2.71




### COSMOSIL Chromatogram Index

Sample: Isopropyl β-D-1-thiogalactopyranoside  
CAS No.: [367-93-1]  
Molecular formula:  $C_{12}H_{21}O_6S$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6, Atten=8  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 6.23min  
Capacity factor: 1.15




### COSMOSIL Chromatogram Index

Sample: Kojic Acid  
CAS No.: [501-30-4]  
Molecular formula:  $C_6H_6O_4$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV245 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.1mg/ml  
Injection volume: 1.0µl  
Retention time: 3.83min  
Capacity factor: 0.46




### COSMOSIL Chromatogram Index

Sample: Malonic Acid  
CAS No.: [141-82-2]  
Molecular formula:  $C_3H_4O_4$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 10.78min  
Capacity factor: 2.81




### COSMOSIL Chromatogram Index

Sample: Mecobalamin  
CAS No.: [13422-55-4]  
Molecular formula:  $C_{50}H_{84}CoN_{14}O_{14}P$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV266 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0µl  
Retention time: 6.22min  
Capacity factor: 1.35




### COSMOSIL Chromatogram Index

Sample: L-Leucine  
CAS No.: [61-90-5]  
Molecular formula:  $C_6H_{13}NO_2$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 4.91min  
Capacity factor: 0.87




### COSMOSIL Chromatogram Index

Sample: D-Leucyl-L-tyrosine  
CAS No.: [3303-29-5]  
Molecular formula:  $C_{13}H_{21}N_3O_4$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 254nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 7.79min  
Capacity factor: 1.96




### COSMOSIL Chromatogram Index

Sample: Metanilic Acid  
CAS No.: [121-47-1]  
Molecular formula:  $C_6H_7NO_3S$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 13.68min  
Capacity factor: 3.80




### COSMOSIL Chromatogram Index

Sample: L-Methionine  
CAS No.: [63-68-3]  
Molecular formula:  $C_5H_{11}NO_2S$   
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)=60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 4.15min  
Capacity factor: 0.54



### COSMOSIL Chromatogram Index


Sample: *N*-Methylglucamine  
CAS No.: [6284-40-8]  
Molecular formula: C<sub>7</sub>H<sub>17</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI  
Attenuation: 4 × 10<sup>9</sup>RIU/FS  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0μl  
Retention time: 15.52min  
Capacity factor: 4.22



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *N*-Methylhydroxylamine  
CAS No.: [593-77-1]  
Molecular formula: CH<sub>3</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELS  
Gain=6, Atten=8  
Sample conc.: 1.0mg/ml  
Injection volume: 2.0μl  
Retention time: 13.45min  
Capacity factor: 4.21



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *DL*-Norleucine  
CAS No.: [616-06-8]  
Molecular formula: C<sub>6</sub>H<sub>13</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0μl  
Retention time: 4.89min  
Capacity factor: 0.86



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *DL*-Norvaline  
CAS No.: [760-78-1]  
Molecular formula: C<sub>6</sub>H<sub>11</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5μl  
Retention time: 5.43min  
Capacity factor: 1.07



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: 6-Methyl-2-thiouracil  
CAS No.: [56-04-2]  
Molecular formula: C<sub>5</sub>H<sub>4</sub>N<sub>2</sub>OS  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5μl  
Retention time: 5.58min  
Capacity factor: 0.84



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Mucic Acid  
CAS No.: [526-99-8]  
Molecular formula: C<sub>6</sub>H<sub>8</sub>O<sub>6</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0μl  
Retention time: 10.27min  
Capacity factor: 2.62



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *L*-Ornithine  
CAS No.: [70-26-8]  
Molecular formula: C<sub>6</sub>H<sub>12</sub>N<sub>2</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)=60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0μl  
Retention time: 8.39min  
Capacity factor: 2.10



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Orotic Acid  
CAS No.: [65-86-1]  
Molecular formula: C<sub>5</sub>H<sub>4</sub>N<sub>2</sub>O<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.5mg/ml  
Injection volume: 1.0μl  
Retention time: 15.24min  
Capacity factor: 4.36



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Murexide  
CAS No.: [3051-09-0]  
Molecular formula: C<sub>8</sub>H<sub>6</sub>N<sub>2</sub>O<sub>6</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5μl  
Retention time: 13.47min  
Capacity factor: 3.69



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Nicotinamide  
CAS No.: [98-92-0]  
Molecular formula: C<sub>6</sub>H<sub>6</sub>N<sub>2</sub>O  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.2mg/ml  
Injection volume: 1.0μl  
Retention time: 5.40min  
Capacity factor: 0.77



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Oxalic Acid  
CAS No.: [144-62-7]  
Molecular formula: C<sub>2</sub>H<sub>2</sub>O<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 5.0mg/ml  
Injection volume: 0.5μl  
Retention time: 12.08min  
Capacity factor: 3.27



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Oxamic Acid  
CAS No.: [471-47-6]  
Molecular formula: C<sub>3</sub>H<sub>4</sub>N<sub>2</sub>O<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.1mg/ml  
Injection volume: 1.0μl  
Retention time: 4.83min  
Capacity factor: 0.71



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Nicotinic Acid  
CAS No.: [59-67-6]  
Molecular formula: C<sub>6</sub>H<sub>5</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5μl  
Retention time: 16.67min  
Capacity factor: 4.87



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *L*-Noradrenaline  
CAS No.: [51-41-2]  
Molecular formula: C<sub>8</sub>H<sub>9</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Citrate buffer(pH7.0)=60/40  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 auFS  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0μl  
Retention time: 5.47min  
Capacity factor: 1.07



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: Oxytocin  
CAS No.: [50-56-6]  
Molecular formula: C<sub>43</sub>H<sub>64</sub>N<sub>12</sub>O<sub>12</sub>S<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 0.4mg/ml  
Injection volume: 0.5μl  
Retention time: 3.71min  
Capacity factor: 0.39



NACALAI TESQUE, INC

### COSMOSIL Chromatogram Index


Sample: *D*-Pantoic Acid  
CAS No.: [79-83-4]  
Molecular formula: C<sub>8</sub>H<sub>15</sub>NO<sub>4</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 auFS  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0μl  
Retention time: 13.21min  
Capacity factor: 3.60



NACALAI TESQUE, INC


### COSMOSIL Chromatogram Index

Sample: *L*-(*-*)-Phenylalanine  
CAS No.: [63-91-2]  
Molecular formula:  $C_9H_9NO_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 254nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 5.10min  
Capacity factor: 0.94




### COSMOSIL Chromatogram Index

Sample: *p*-Phenylenediamine  
CAS No.: [106-50-3]  
Molecular formula:  $C_6H_8N_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 4.15min  
Capacity factor: 0.36




### COSMOSIL Chromatogram Index

Sample: *L*-Proline  
CAS No.: [147-85-3]  
Molecular formula:  $C_5H_9NO_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.35min  
Capacity factor: 1.22




### COSMOSIL Chromatogram Index

Sample: Propionic Acid  
CAS No.: [79-09-4]  
Molecular formula:  $C_3H_6O_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 14.85min  
Capacity factor: 4.24




### COSMOSIL Chromatogram Index

Sample: *L*-(+)- $\alpha$ -Phenylglycine  
CAS No.: [2935-35-5]  
Molecular formula:  $C_9H_9NO_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0µl  
Retention time: 5.96min  
Capacity factor: 1.27




### COSMOSIL Chromatogram Index

Sample: Phosphocreatine  
CAS No.: [67-07-2]  
Molecular formula:  $C_4H_9N_3O_5P$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 11.42min  
Capacity factor: 3.00




### COSMOSIL Chromatogram Index

Sample: Pyridoxine  
CAS No.: [65-23-6]  
Molecular formula:  $C_8H_{10}NO_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5µl  
Retention time: 6.78min  
Capacity factor: 1.35




### COSMOSIL Chromatogram Index

Sample: Pyruvic Acid  
CAS No.: [127-17-3]  
Molecular formula:  $C_3H_4O_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 2.0µl  
Retention time: 3.97min  
Capacity factor: 0.39




### COSMOSIL Chromatogram Index

Sample: *O*-Phospho-*L*-serine  
CAS No.: [407-41-0]  
Molecular formula:  $C_4H_8NO_6P$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 3.0µl  
Retention time: 9.19min  
Capacity factor: 2.24




### COSMOSIL Chromatogram Index

Sample: Picolinic acid  
CAS No.: [98-98-6]  
Molecular formula:  $C_5H_5NO_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Citrate buffer(pH7.0)=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 5.03min  
Capacity factor: 0.92




### COSMOSIL Chromatogram Index

Sample: Ribose-5-phosphate  
CAS No.: [4306-28-1]  
Molecular formula:  $C_5H_{10}O_8P$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: RI  
Attenuation:  $4 \times 10^3$  RIU/FS  
Sample conc.: 10.0mg/ml  
Injection volume: 5.0µl  
Retention time: 8.02min  
Capacity factor: 2.06




### COSMOSIL Chromatogram Index

Sample: *D*-Saccharic Acid  
CAS No.: [87-73-0]  
Molecular formula:  $C_6H_{12}O_8$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 10.48min  
Capacity factor: 2.69




### COSMOSIL Chromatogram Index

Sample: Pivalic Acid  
CAS No.: [75-98-9]  
Molecular formula:  $C_5H_{10}O_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 12.14min  
Capacity factor: 3.28




### COSMOSIL Chromatogram Index

Sample: Procaterol  
CAS No.: [72332-33-3]  
Molecular formula:  $C_{16}H_{23}N_3O_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 20mmol/l Ammonium acetate =85/15  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 6.17min  
Capacity factor: 1.25




### COSMOSIL Chromatogram Index

Sample: Sarcosine  
CAS No.: [107-97-1]  
Molecular formula:  $C_4H_9NO_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 210nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0µl  
Retention time: 6.30min  
Capacity factor: 1.40




### COSMOSIL Chromatogram Index

Sample: Sebacic Acid  
CAS No.: [111-20-6]  
Molecular formula:  $C_{18}H_{34}O_4$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 1.5µl  
Retention time: 9.43min  
Capacity factor: 2.28




### COSMOSIL Chromatogram Index

Sample: L-Serine  
CAS No.: [56-45-1]  
Molecular formula:  $C_3H_7NO_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 7.92min  
Capacity factor: 2.01




### COSMOSIL Chromatogram Index

Sample: Sinigrin  
CAS No.: [3952-98-5]  
Molecular formula:  $C_{16}H_{19}KNO_8S_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 12.57min  
Capacity factor: 3.38




### COSMOSIL Chromatogram Index

Sample: L-Theanine  
CAS No.: [3081-61-6]  
Molecular formula:  $C_7H_{14}N_2O_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 0.5µl  
Retention time: 5.89min  
Capacity factor: 1.21




### COSMOSIL Chromatogram Index

Sample: 2-Thiobarbituric Acid  
CAS No.: [504-17-6]  
Molecular formula:  $C_4H_4N_2O_3S$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 5.18min  
Capacity factor: 0.82




### COSMOSIL Chromatogram Index

Sample: Sorbic Acid  
CAS No.: [110-44-1]  
Molecular formula:  $C_8H_8O_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 13.19min  
Capacity factor: 3.59




### COSMOSIL Chromatogram Index

Sample: Succinic Acid  
CAS No.: [110-15-6]  
Molecular formula:  $C_4H_6O_4$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 0.5µl  
Retention time: 10.64min  
Capacity factor: 2.74




### COSMOSIL Chromatogram Index

Sample: 2-Thiouacil  
CAS No.: [141-90-2]  
Molecular formula:  $C_5H_6N_2OS$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5µl  
Retention time: 6.38min  
Capacity factor: 1.11




### COSMOSIL Chromatogram Index

Sample: L-Threonine  
CAS No.: [72-19-5]  
Molecular formula:  $C_4H_9NO_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0µl  
Retention time: 7.19min  
Capacity factor: 1.73




### COSMOSIL Chromatogram Index

Sample: Sulbactam  
CAS No.: [68373-14-8]  
Molecular formula:  $C_8H_{11}NO_5S$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 0.5µl  
Retention time: 10.86min  
Capacity factor: 2.81




### COSMOSIL Chromatogram Index

Sample: Sulfamic acid  
CAS No.: [121-57-3]  
Molecular formula:  $C_6H_7NO_3S$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.10mg/ml  
Injection volume: 1.0µl  
Retention time: 13.87min  
Capacity factor: 3.87




### COSMOSIL Chromatogram Index

Sample: Todralazine  
CAS No.: [14679-73-3]  
Molecular formula:  $C_{11}H_{14}N_2O_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV240 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 4.56min  
Capacity factor: 0.51




### COSMOSIL Chromatogram Index

Sample: Trichloroacetic Acid  
CAS No.: [76-03-9]  
Molecular formula:  $C_2HCl_3O_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 11.83min  
Capacity factor: 3.17




### COSMOSIL Chromatogram Index

Sample: L-(+)-Tartaric Acid  
CAS No.: [87-69-4]  
Molecular formula:  $C_4H_6O_6$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.5µl  
Retention time: 10.52min  
Capacity factor: 2.70




### COSMOSIL Chromatogram Index

Sample: Taurine  
CAS No.: [107-35-7]  
Molecular formula:  $C_2H_7NO_3S$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6,Atten=8  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0µl  
Retention time: 6.40min  
Capacity factor: 1.25




### COSMOSIL Chromatogram Index

Sample: Tris(hydroxymethyl)aminomethane  
CAS No.: [77-86-1]  
Molecular formula:  $C_4H_{11}NO_3$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =80/20  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: ELSD  
Attenuation: Gain=6,Atten=8  
Sample conc.: 2.0mg/ml  
Injection volume: 1.0µl  
Retention time: 6.47min  
Capacity factor: 1.48




### COSMOSIL Chromatogram Index

Sample: L-Tryptophan  
CAS No.: [73-22-3]  
Molecular formula:  $C_{11}H_{12}N_2O_2$   
Column: HILIC  
Column size: 4.6mm.I.D.-250mm  
Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH4.0) =70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.5mg/ml  
Injection volume: 0.5µl  
Retention time: 5.14min  
Capacity factor: 0.95








コスモシール使用者様の声 HILIC 抜粋

"こんなに役立つ情報"


**COSMOSIL Chromatogram Index**

Sample: L-Tyrosine  
CAS No.: [60-18-4]  
Molecular formula: C<sub>9</sub>H<sub>9</sub>NO<sub>3</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV 254nm  
Attenuation: 0.128 au/s  
Sample conc.: 5.0mg/ml  
Injection volume: 1.0μl  
Retention time: 5.92min  
Capacity factor: 1.25




**COSMOSIL Chromatogram Index**

Sample: Uracil  
CAS No.: [66-22-8]  
Molecular formula: C<sub>4</sub>H<sub>4</sub>N<sub>2</sub>O<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 0.5μl  
Retention time: 5.33min  
Capacity factor: 0.84




**COSMOSIL Chromatogram Index**

Sample: Urea  
CAS No.: [57-13-6]  
Molecular formula: CH<sub>4</sub>N<sub>2</sub>O  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 2.0μl  
Retention time: 6.12min  
Capacity factor: 1.15




**COSMOSIL Chromatogram Index**

Sample: Uridine  
CAS No.: [58-96-8]  
Molecular formula: C<sub>9</sub>H<sub>9</sub>N<sub>3</sub>O<sub>6</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV260 nm  
Attenuation: 0.128 au/s  
Sample conc.: 0.1mg/ml  
Injection volume: 1.0μl  
Retention time: 8.30min  
Capacity factor: 1.86




**COSMOSIL Chromatogram Index**

Sample: L-Valine  
CAS No.: [72-18-4]  
Molecular formula: C<sub>6</sub>H<sub>11</sub>NO<sub>2</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV210 nm  
Attenuation: 0.128 au/s  
Sample conc.: 10.0mg/ml  
Injection volume: 1.0μl  
Retention time: 5.63min  
Capacity factor: 1.14




**COSMOSIL Chromatogram Index**

Sample: Vitamin B<sub>1</sub> (Thiamine)  
CAS No.: [67-03-8]  
Molecular formula: C<sub>12</sub>H<sub>17</sub>ClN<sub>4</sub>OS  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/50mmol/l Ammonium acetate=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 1.0μl  
Retention time: 11.25min  
Capacity factor: 2.93



**COSMOSIL Chromatogram Index**

Sample: Vitamin B<sub>2</sub> (Riboflavin)  
CAS No.: [83-88-5]  
Molecular formula: C<sub>17</sub>H<sub>20</sub>N<sub>4</sub>O<sub>6</sub>  
Column: HILIC  
Column size: 4.6mm I.D.-250mm  
Mobile phase: Acetonitrile/10mmol/l Ammonium acetate=90/10  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220 nm  
Attenuation: 0.128 au/s  
Sample conc.: 1.0mg/ml  
Injection volume: 0.5μl  
Retention time: 11.33min  
Capacity factor: 2.92



ナカライテスクでは、C<sub>18</sub> カラムの分離で困った時のカラムとしてコスモシール高機能カラム Cholester、 NAP、 HILIC、 PYE、 PBB-R を紹介しています。  
今回、コスモシールを使用されていたお客様からいただきました "こんなに役立つ情報" から HILIC についてのお客様の声ををご紹介します。(研究者の皆様のご所属等の情報は、データご提供時の情報に基づいています。)

HILIC 情報抜粋

提供者名	近畿大学薬学総合研究所 食品薬学研究室 森川 敏生 准教授
カラムサイズ	4.6 × 250, 20 × 250
サンプル名	天然由来化合物 (配糖体成分)
内容	今まで C <sub>18</sub> カラムで分離が出来なかったアシル基の結合位置の異なる配糖体成分がコスモシール HILIC を使用することにより完全分離する事が出来た。
提供者名	大学 薬学部 研究者様
カラムサイズ	4.6 × 250, 10 × 250
サンプル名	糖の結合する位置が異なるフラボノイド C 配糖体の分離
内容	結合する位置が異なるフラボノイド C 配糖体の分離を試みたとき、ODS はだめでゲル濾過カラム等ではやはり分子量も同じなので分離はできなかったのですが、コスモシール HILIC を使ったら一発でした。(条件: 98% CH <sub>3</sub> CN)
提供者名	食品メーカー 研究者様
カラムサイズ	2.0 × 150
サンプル名	脂質、糖質、ステロイド類、脂溶性ビタミンなど
内容	コスモシール HILIC は他社のカラムと比べ、分離が良く使いやすい。MS 用のカラムとして十分利用できる。
提供者名	大学 薬学部 研究者様
カラムサイズ	4.6 × 250
サンプル名	微生物 Pseudomonas 属が分泌する細胞外代謝産物
内容	非常に親水性が高い化合物であるため C <sub>18</sub> カラムによる分析が不可能でした。コスモシール HILIC を用いることで、再現性のある分析が行うことができました。
提供者名	公的研究機関 研究者
カラムサイズ	2.0 × 50, 4.6 × 250, 10 × 250
サンプル名	化粧品など
内容	コスモシール HILIC を使用することで既報の方法ではうまく HPLC 分離できなかった成分を分離できるようになりました。さらには、試料中で生じていた分解物を分取精製し、構造決定にまで至るなど、非常に重宝しています。今後の課題として、耐久性の向上に期待しています。

"こんなに役立つ情報" 提供依頼

今回ご紹介しましたようなコスモシール使用者様の声を "こんなに役立つ情報" として充実・共有化することで、皆様の研究にお役に立ちたいと考えています。  
今後、"こんなに役立つ情報" を充実させることで、使用者様の声が共有化できるコミュニティサイトのオープンを検討したいと考えていますので、是非とも、使用者様の "こんなに役立つ情報" をお寄せください。

(投稿アドレス)

弊社 Web に投稿ページを準備していますので、是非、投稿を宜しくお願い申し上げます。  
<http://www.nacalai.co.jp/update/cosmosil.html>

Technical Note

データご提供：近畿大学 薬学総合研究所 食品薬学研究室 森川 敏生 准教授


クロマトグラフィーで汎用されている C<sub>18</sub> カラムは、疎水性に差がないような糖置換体や糖の結合位置が異なる化合物に対して分離が不十分な場合があります。このような化合物の中で親水性基に差があるような化合物に対して親水性相互作用を有するコスモシール HILIC を用いることにより分離が改善する場合があります。今回、アシル基の結合位置が異なる配糖体の分析データをご提供いただきましたのでご紹介いたします。

(1) 結論

キク科 (Asteraceae) 植物ダイジー (Bellis perennis) は西ヨーロッパ原産の多年生草本で、その花期が長いことから和名をヒナギク (雛菊)、チョウメイギク (長命菊) およびエンメイギク (延命菊) などと称されます。日本ではおもに園芸品種として親しまれていますが、西欧諸国においては、古くからその開花時の全草や根部を打撲傷、出血、筋肉痛および皮膚病やリウマチの治療に用いられているとともに、その若葉や蕾、花弁などをサラダに加えて生食するなど、食用としても供されています。また、ダイジーの花部であるダイジーフラワーはハーブティーなどにも利用され、高度に酸化されたアグリコン部 (非糖部) を有するビスデスモシド型トリテルペンサポニン類が含有成分として報告されています。とりわけ、アシル基の有無や種類、あるいは結合位置の違いなど、その化学構造が微妙に異なるサポニン類が多数報告されています。

(2) アシル基の結合位置の異なるサポニン成分の分析


一般的に親水性相互作用クロマトグラフィー用カラムは、疎水性が非常に小さく C<sub>18</sub> カラムで保持が小さいような化合物に対して有効なカラムとなります。C<sub>18</sub> カラムで十分保持があるアシル基の結合位置の異なる化合物に対して C<sub>18</sub> カラムおよびコスモシール HILIC を用いて分析を行いました。



COSMOSIL Application Data



Column: C<sub>18</sub> column (Compeitor)  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ Methanol/ H<sub>2</sub>O = 35/16/49 + 1%Acetic Acid  
 Flow rate: 1.0 ml/min  
 Temperature: room temperature  
 Detection: RI

Component: 1; Perennisaponin J (1)  
 2; Perennisaponin K (2)



Column: COSMOSIL HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 1%Acetic Acid- H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: room temperature  
 Detection: UV220nm

Component: 1; Perennisaponin J (1)  
 2; Perennisaponin K (2)



Data courtesy of Toshio Morikawa, Ph.D. Associate Professor, Pharmaceutical Food Sciences, Pharmaceutical Research and Technology Institute, Kinki University

NACALAI TESQUE, INC

結果、疎水性が類似しているため C<sub>18</sub> カラムで分離できませんでしたが、コスモシール HILIC を用いることによりアシル基の結合位置のみが異なるサポニン (perennisaponin J および K) の分離について、極めて良好な分離ができました。


(3) コスモシール HILIC による分取精製

内径 20mm の分取カラムを用いて、それぞれのサポニンを精製しました。

COSMOSIL Application Data

Column: COSMOSIL HILIC  
 Column size: 20mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 1%Acetic Acid-H<sub>2</sub>O = 90/10  
 Flow rate: 0.9 ml/min  
 Temperature: room temperature  
 Detection: UV220nm

Component: 1; Perennisaponin J (1)  
 2; Perennisaponin K (2)



Data courtesy of Toshio Morikawa, Ph.D. Associate Professor, Pharmaceutical Food Sciences, Pharmaceutical Research and Technology Institute, Kinki University

結果、コスモシール HILIC においてセミ分取サイズのカラムにおいても再現し、それぞれのサポニンが単離することができました。

(4) 考察

C<sub>18</sub> カラムで分離できなかったアシル基の結合位置のみが異なるサポニン (perennisaponin J および K) の分離について、コスモシール HILIC を用いて検討したところ極めて良好な分離ができました。この良好な分離はセミ分取サイズのカラムにおいても再現し、それぞれのサポニンが単離することができ、NMR などによる構造解析および各種フィジカルデータ測定に必要なサンプル量を確保することができました。

(5) 参考文献

T.Morikawa et.al Acylated Oleanane-Type Triterpene Bisdesmosides: Perennisaponins G, H, I, J, K, L, and M with Pancreatic Lipase Inhibitory Activity from the Flowers of Bellis perennis. Helv. Chim. Acta 93 573-586 (2010).

年々親水性相互作用クロマトグラフィー関連の文献は増加傾向にあり、注目されている手法の一つです。ここでは、弊社オリジナルカラムであるコスモシルHILICが記載されている最近の文献リストの一部をご紹介します。なお、著作権の関係で弊社より文献をご提供することはできませんのでご了承ください。

## HILIC 文献抜粋

No.	Title	AUTHOR	JOURNAL	ISSUE	PAGE	YEAR
1	Characterization of the decomposition of compounds derived from imidazolidinyl urea in cosmetics and patch test materials	Takahiro Doi, Akihiro Takeda, Akiko Asada, Keiji Kajimura	Contact Dermatitis	67 (5)	284-292	2012
2	A simple graphical representation of selectivity in hydrophilic interaction liquid chromatography	Mohammed E.A. Ibrahim, Yang Liu, Charles A. Lucy	Journal of Chromatography A	1260	126-131	2012
3	Comparison of 2-amino-[3- <sup>11</sup> C] isobutyric acid and 2-deoxy-2-[ <sup>18</sup> F] fluoro-D-glucose in nude mice with xenografted tumors and acute inflammation	Tsuiji, Atsushi B; Kato, Koichi; Sugyo, Aya; Okada, Makki; Sudo, Hitomi; Yoshida, Chisato; Wakizaka, Hidekatsu; Zhang, Ming-Rong; Saga, Tsuneo	Nuclear Medicine Communications	33 (10)	1058-1064	2012
4	In Vitro and in Vivo Metabolism of Verproside in Rats	Min Gi Kim, Deok-Kyu Hwang, Hyeon-Uk Jeong, Hye Young Ji, Sei-Ryang Oh, Yongnam Lee, Ji Seok Yoo, Dae Hee Shin and Hye Suk Lee	Molecules	17 (10)	11990-12002	2012
5	Functional expression of carnitine/organic cation transporter OCTN1 in mouse brain neurons: Possible involvement in neuronal differentiation	Noritaka Nakamichi, Takayuki Taguchi, Hiroshi Hosotani, Tomohiko Wakayama, Takuya Shimizu, Tomoko Sugiura, Shoichi Iseki, Yukio Kato,	Neurochemistry International	In Press		2012
6	Characterization and use of hydrophilic interaction liquid chromatography type stationary phases in supercritical fluid chromatography	Caroline West, Syame Khater, Eric Lesellier	Journal of Chromatography A	1250	182-195	2012
7	Analysis of 8-hydroxy-2-deoxyguanosine in human urine using hydrophilic interaction chromatography with tandem mass spectrometry	Chiemi Hosozumi, Akira Toriba, Thanyarat Chuesaard, Takayuki Kameda, Ning Tang, Kazuichi Hayakawa	Journal of Chromatography B	893-894	173-176	2012
8	A NOVEL NORMAL PHASE HPLC METHOD FOR THE QUANTIFICATION OF N-FORMYL IMPURITY IN AZACITIDINE ACTIVE PHARMACEUTICAL INGREDIENTS AND PHARMACEUTICAL DOSAGE FORMS	T. Satyanarayana Raju, L. Kalyanaraman, K. S. V. Raghavachary & P. Yadagiri Swamy	Journal of Liquid Chromatography & Related Technologies	35 (8)	1070-1080	2012
9	Determination of Histamine in Seafood by Hydrophilic Interaction Chromatography/Tandem Mass Spectrometry	Tatsuo YOSHIDA, Hirotohi HAMADA, Hiroshi MURAKAWA, Hidekazu YOSHIMOTO, Toshiaki TOBINO, Kei TODA	Analytical Sciences	28 (2)	179-182	2012
10	Triazole-Linked DNA as a Primer Surrogate in the Synthesis of First-Strand cDNA	Dr. Tomoko Fujino, Dr. Ken-ichi Yasumoto, Naomi Yamazaki, Ai Hasome, Prof. Kazuhiro Sogawa, Prof. Hiroyuki Isebe	Chemistry An Asian Journal	6 (11)	2956-2960	2011
11	Retention and selectivity of stationary phases for hydrophilic interaction chromatography	Yong Guo, Sheetal Gaiki	Journal of Chromatography A	1218 (35)	5920-5938	2011
12	Chromatographic characterization of hydrophilic interaction liquid chromatography stationary phases: Hydrophilicity, charge effects, structural selectivity, and separation efficiency	Yuusuke Kawachi, Tohru Ikegami, Hirota Takubo, Yuka Ikegami, Masatoshi Miyamoto, Nobuo Tanaka	Journal of Chromatography A	1218 (35)	5903-5919	2011
13	The different decomposition properties of diazolidinyl urea in cosmetics and patch test materials	Takahiro Doi, Keiji Kajimura, Shuzo Taguchi	Contact Dermatitis	65 (2)	81-91	2011
14	Stationary and mobile phases in hydrophilic interaction chromatography: a review	Pavel Jandera	Analytica Chimica Acta	692 (1-2)	1-25	2011
15	Degradation of N-Acetyl-D-glucosamine and D-Glucosamine in Subcritical Water and Properties of the Degradation Products	Rongchun WANG, Takashi KOBAYASHI and Shuji ADACHI	Food Science and Technology Research	17 (4)	273-278	2011
16	Determination of isoascorbic acid in fish tissue by hydrophilic interaction liquid chromatography ultraviolet detection	Spyros Drivelos, Marilena E. Dasenaki and Nikolaos S. Thomaidis	Analytical and Bioanalytical Chemistry	397 (6)	2199-2210	2010

No.	Title	AUTHOR	JOURNAL	ISSUE	PAGE	YEAR
17	Hepatoprotective Effects of Flavonoids from Shekwasha (Citrus depressa) against D-Galactosamine-Induced Liver Injury in Rats	Toshiyuki AKACHI, Yasuyuki SHIINA, Yaiyo OHISHI, Takumi KAWAGUCHI, Hirokazu KAWAGISHI, Tatsuya MORITA, Makoto MORI and Kimio SUGIYAMA	J. Nutr Sci Vitaminol	56 (1)	60-67	2010
18	A Novel Glucosylation Reaction on Anthocyanins Catalyzed by Acyl-Glucose Dependent Glucosyltransferase in the Petals of Carnation and Delphinium	Yuki Matsuba, Nobuhiro Sasaki, Masayuki Tera, Masachika Okamura, Yutaka Abe, Emi Okamoto, Haruka Nakamura, Hisakage Funabashi, Makoto Takatsu, Mikako Saito, Hideaki Matsuoka, Kazuo Nagasawa and Yoshihiro Ozekia	The Plant Cell	22 (10)	3374-3389	2010
19	Molecular identification of unsaturated uronate reductase prerequisite for alginate metabolism in <i>Spingomonas</i> sp. A1	Ryuichi Takase, Akihito Ochiai, Bunzo Mikami, Wataru Hashimoto, Kousaku Murata,	Biochimica et Biophysica Acta (BBA) - Proteins & Proteomics	1804 (9)	1925-1936	2010
20	Inhibitory Effects of Acylated Acyclic Sesquiterpene Oligoglycosides from the Pericarps of <i>Sapindus rarak</i> on Tumor Necrosis Factor-Induced Cytotoxicity	Toshio Morikawa, Yuanyuan Xie, Kiyofumi Ninomiya, Masaki Okamoto, Osamu Muraoka, Dan Yuan, Masayuki Yoshikawa and Takao Hayakawa	Chem. Pharm. Bull.	58 (9)	1276-1280	2010
21	Approach to hydrophilic interaction chromatography column selection: Application to neurotransmitters analysis	Raluca-Ioana Chirita, Caroline West, Adriana-Luminita Finaru, Claire Elfakir	Journal of Chromatography A	1217 (18)	3091-3104	2010
22	Medicinal Flowers. Part 29: Acylated Oleanane-Type Triterpene Bidesmosides: Perennisaponins G, H, I, J, K, L, and M with Pancreatic Lipase Inhibitory Activity from the Flowers of <i>Bellis perennis</i>	Toshio Morikawa, Xuezheng Li, Eriko Nishida, Seikou Nakamura, Kiyofumi Ninomiya, Hisashi Matsuda, Yoshimi Oda, Osamu Muraoka, Masayuki Yoshikawa	Helvetica Chimica Acta	93 (3)	573-586	2010
23	Unusual amino acid derivatives from the mushroom <i>Pleurocybella porrigens</i>	Takumi Kawaguchi, Tomohiro Suzuki, Yuka Kobayashi, Shinya Kodani, Hirofumi Hirai, Kaoru Nagai, Hirokazu Kawagishi	Tetrahedron	66 (2)	504-507	2010
24	Structures of Acetylated Oleanane-Type Triterpene Saponins, Rarasaponins IV, V, and VI, and Anti-hyperlipidemic Constituents from the Pericarps of <i>Sapindus rarak</i>	Yasunobu Asao, Toshio Morikawa, Yuanyuan Xie, Masaki Okamoto, Makoto Hamao, Hisashi Matsuda, Osamu Muraoka, Dan Yuan and Masayuki Yoshikawa	Chem. Pharm. Bull.	57 (2)	198-203	2009
25	Development and validation of a reversed-phase high-performance liquid chromatographic method for quantification of peptide dendrimers in human skin permeation experiments	S. Mutalik, A.K. Hewavitharana, P.N. Shaw, Y.G. Anissimov, M.S. Roberts, H.S. Parekh,	Journal of Chromatography B	877 (29)	3556-3562	2009
26	Determination of para-aminohippuric acid (PAH) in human plasma and urine by liquid chromatography tandem mass spectrometry	Phey Yen Han, P. Nicholas Shaw, Carl M.J. Kirkpatrick	Journal of Chromatography B	877 (27)	3215-3220	2009
27	Oxidation of Methionine to Dehydromethionine by Reactive Halogen Species Generated by Neutrophils	Alexander V. Peskin, Rufus Turner, Ghassan J. Maghzal, Christine C. Winterbourn and Anthony J. Kettle	Biochemistry	48 (42)	10175-10182	2009
28	Oleanane-type triterpene oligoglycosides with pancreatic lipase inhibitory activity from the pericarps of <i>Sapindus rarak</i>	Toshio Morikawa, Yuanyuan Xie, Yasunobu Asao, Masaki Okamoto, Chihiro Yamashita, Osamu Muraoka, Hisashi Matsuda, Yutana Pongpiriyadacha, Dan Yuan, Masayuki Yoshikawa	Phytochemistry	70 (9)	1166-1172	2009
29	Direct Evidence for Efficient Transport and Minimal Metabolism of L-Cephalexin by Oligopeptide Transporter 1 in Budded Baculovirus Fraction	Keisuke Mitsuoka, Ikumi Tamai, Yasushi Morohashi, Yoshiyuki Kubo, Ryoichi Saitoh, Akira Tsuji and Yukio Kato	Biol. Pharm. Bull.	32 (8)	1459-1461	2009
30	Simultaneous measurement of diazolidinyl urea, urea, and allantoin in cosmetic samples by hydrophilic interaction chromatography	Takahiro Doi, Keiji Kajimura, Satoshi Takatori, Naoki Fukui, Shuzo Taguchi, Shozo Iwagami	Journal of Chromatography B	877 (10)	1005-1010	2009
31	親水性相互作用クロマトグラフィー (HILIC) 開発の背景、および分離モードの特徴	池上亨, 田窪宏貴, 田中信男	Chromatography	29 (2)		2008
32	Tetrodotoxin poisoning evidenced by solid-phase extraction combining with liquid chromatography tandem mass spectrometry	Hsiao-Chin Jen, Shin-Jung Lin, Yung-Hsiang Tsai, Chun-Hsiang Chen, Zu-Chun Lin, Deng-Fwu Hwang,	Journal of Chromatography B	871 (1)	95-100	2008

	Sample name	Page	
A	Acesulfame	10	
	Acetamide	12	
	Acetazolamide	12	
	Acetizonic Acid	12	
	Acrylic Acid	12	
	L- Alanine	12	
	-Alanine	12	
	Allantoic Acid	11	
	Allantoin	11, 12	
	p-Aminobenzamide	12	
	p-Aminobenzoic Acid	13	
	4-Amino-n-butyrilic Acid [ GABA ]	13	
	6-Aminohexanoic Acid [ 6-Amino-n-caproic Acid ]	13	
	5-Aminolevulinic Acid	13	
	2-Aminopyridine	13	
	3-Aminopyridine	13	
	5-Amino-1H-tetrazole	13	
	3-Amino-1H-1,2,4-triazole	13	
	5-Aminouracil	14	
	Ammelide	8	
	Ammeline	8	
	Amphotericin B	14	
	Angiotensin I (Human)	5, 14	
	Angiotensin II (Human)	9	
	Angiotensin II, [ Asn <sup>1</sup> ,Val <sup>6</sup> ]	5, 8, 14	
	Angiotensin II, [ Sar <sup>1</sup> ,Ala <sup>8</sup> ]	5, 14	
	Angiotensin II, [ Sar <sup>1</sup> ,Ile <sup>6</sup> ]	5, 8, 14	
	Angiotensin II, [ Sar <sup>1</sup> ,Thr <sup>6</sup> ]	5, 14	
	Angiotensin II, [ Val <sup>6</sup> ]	5, 8, 15	
	Angiotensin II, Des-Asp <sup>1</sup> - [ Ile <sup>6</sup> ]	5, 8, 15	
	L-Arginine	15	
	Ascorbic Acid	4, 7, 8	
	L-(+)-Ascorbic Acid [ Vitamin C ]	8, 15	
	L-Asparagine	15	
	Aspartame	10	
	L-Aspartic Acid	15	
	6-Azauracil	15	
	Aztreonam	15	
	B	Benzamide	16
		Benzenesulfonic Acid	16
		Benzoic Acid	10, 16
Benzylamine		16	
Bromoacetic Acid	16		
C	Cacotheline	16	
	Ca eine	10	
	Camostat	16	
	L-Carnitine	16	
	Ceftriaxone	17	
	Chloroacetic Acid	17	
	Citrazinic Acid	17	
	L-Citrulline	8	
	Choline Chloride	11	
	C	Choline Hydrogen Tartrate	11
		Creatine	17
		Creatinine	10, 17
		Cyanoacetic Acid	17
		Cyanuric Acid	8, 17
		L-Cysteine	17
		L-(-)-Cystine	18
Cytidine		18	
Cytosine		18	
D		2'-Deoxyguanosine	10
		3,4-Diaminobenzoic Acid	18
		3,5-Diaminobenzoic Acid	18
		2,4-Diaminophenol	18
		DL-2,6-Diaminopimelic Acid	18
		DL-2,3-Diaminopropionic Acid	18
		Diatrizoic Acid	19
		Diethylene Glycol	8
		Dipicolinic Acid	19
		Dithiouracil	19
		L-DOPA	19
Dopamine		19	
E		L-(+)-Ergothioneine	19
		meso-Erythritol	3, 9, 19
Ethylene Glycol	9		
F	Famotidine	19	
	Folic Acid	20	
	Folinic Acid	20	
	Formamide	20	
	D-Fructose-6-phosphate	9, 20	
	Fuchsine, Acid [ Rubin S ]	20	
	Fumaric Acid	20	
	G	GABA [ 4-Amino-n-butyrilic Acid ]	13
		Gluconic Acid	20
		Glucose	9
D-Glucose-6-phosphate		9, 20	
-D-Glucose-1-phosphate		9, 21	
D-Glucuronic Acid		21	
L-Glutamic Acid		21	
L-Glutamine		21	
Glutaric Acid		21	
Glutathione(Reduced Form)		21	
Glyceric Acid		3,9	
DL-Glyceric Acid		21	
Glycerol		8,9	
Glycinamide	21		
Glycine	2, 9, 22		
Glycolic Acid	22		
Glycylglycine	2, 9, 22		
Guanidoacetic Acid	22		
H	1,2,6-Hexanetriol	22	
	L-Histidine	22	
	L-Homocystine	22	

	Sample name	Page	
H	L-Homoserine	22	
	Hydantoic Acid	23	
	Hydantoin	23	
	Hydroxylamine-O-sulfonic Acid	23	
	8-Hydroxy-2'-Deoxyguanosine	10	
	8-Hydroxy Guanosine	10	
	cis-4-Hydroxy-D-proline	23	
	L-Hydroxyproline [ trans-4-Hydroxy-L-proline ]	23	
	N-Hydroxysuccinimide	23	
	I	Indigo carmine	23
		Isoascorbic Acid [ Erythorbic Acid ]	4, 7, 8
		D-Isoascorbic Acid	23
		Isocinchomeronic Acid [ Pyridine-2,5-dicarboxylic Acid ]	24
Isoleucine		9	
L-Isoleucine		24	
Isonicotinic Acid		24	
Isonicotinohydrazide		24	
Isopropyl - -D-1-thiogalactopyranoside [ IPTG ]		24	
Kojic Acid		24	
L	Leucine	9	
	L-Leucine	24	
	D-Leucyl-L-tyrosine	24	
L-Lysine	25		
M	Melamine	4, 8, 25	
	Maleic Acid	25	
	Malic Acid	8	
	L-(-)-Malic Acid	25	
	Malonic Acid	25	
	Mecobalamin	25	
	Metanilic Acid	25	
	L-Methionine	25	
	N-Methylglucamine	26	
	N-Methylhydroxylamine	26	
2-Methylimidazole	10		
4-Methylimidazole	10		
6-Methyl-2-thiouracil	26		
Mucic Acid	26		
Murexide	26		
N	1,5-Naphthalenedisulfonic Acid	11	
	Nicotinamide	8, 26	
	Nicotinic Acid	8, 26	
	L-Noradrenaline [ Norepinephrine ]	26	
	Norepinephrine [ L-Noradrenaline ]	26	
	DL-Norleucine	27	
	DL-Norvaline	27	
	L-Ornithine	27	
O	Orotic Acid	27	
	Oxalic Acid	3, 9, 27	
	Oxamic Acid	3, 9, 27	
	Oxytocin	27	
	P	D-Pantothenic Acid	8, 27
Perennisaponin J		34, 35	
P	Perennisaponin K	34, 35	
	L-(-)-Phenylalanine	28	
	p-Phenylenediamine	28	
	L-(+)- -Phenylglycine	28	
	Phosphocreatine	28	
	O-Phospho-L-serine	28	
	Picolinic Acid	28	
	Pivalic Acid	28	
	Procaterol	28	
	L-Proline	29	
	Propionic Acid	29	
Pyridoxine [ Vitamin B <sub>6</sub> ]	8, 29		
Pyruvic Acid	29		
Q	Quinine	10	
R	Riboflavin [ Vitamin B <sub>2</sub> ]	8, 32	
	Ribose-5-phosphate	29	
S	D-Saccharic Acid	29	
	Saccharin	10	
	Sarcosine	29	
	Sebacic Acid	29	
	L-Serine	30	
	Sinigrin	30	
	Sorbic Acid	8, 10, 30	
	Succinic Acid	30	
	Sulbactam	30	
	Sulfanilic Acid	30	
T	L-(-)-Tartaric Acid	30	
	Taurine	10, 30	
	L-Theanine	31	
	Thiamine [ Vitamin B <sub>1</sub> ]	32	
	2-Thiobarbituric Acid	31	
	2-Thiouracil	31	
	L-Threonine	31	
	Todralazine	31	
	Trichloroacetic Acid	31	
	Trimethylene Glycol	9	
	Tris(hydroxymethyl)aminomethane	3, 9, 31	
Tryptophan	7		
L-Tryptophan	31		
L-Tyrosine	32		
U	Uracil	10, 32	
	Urea	10, 32	
Uridine	10, 32		
V	Valine	9	
	L-Valine	28	
	Vitamin B <sub>1</sub> [ Thiamine ]	28	
	Vitamin B <sub>2</sub> [ Riboflavin ]	8, 32	
	Vitamin B <sub>6</sub> [ Pyridoxine ]	8, 29	
	Vitamin C [ L-(+)-Ascorbic Acid ]	8, 15	