INVIVOGEN'S PRODUCTS FOCUS

LPS

HIGH-QUALITY LPS WITH TWO GRADES OF PURITY

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<https://www.invivogen.com/tlr4-ligands>

“TLR4-specific” LPS from Brand E is contaminated by TLR2 response


Fig.1 HEK-Blue™ hTLR2 cells were stimulated with LPS-EB and LPS-B5 from Brand E, LPS-EB and LPS-B5 Ultrapure, and Pam3CSK4 from InvivoGen at different concentrations. After incubation, NF- κ B-induced SEAP activity was determined by reading the OD value at 655 nm.

TLR4 potency is comparative to the other brand


Fig.2 HEK-Blue™ hTLR4 cells were stimulated with LPS-EB from Brand E and LPS-EB Ultrapure from InvivoGen at different concentrations. After incubation, NF- κ B-induced SEAP activity was determined by reading the OD value at 655 nm.


Fig.3 HEK-Blue™ hTLR4 cells were stimulated with LPS-B5 from Brand E and LPS-B5 Ultrapure from InvivoGen at different concentrations. After incubation, NF- κ B-induced SEAP activity was determined by reading the OD value at 655 nm.

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SMOOTH (S)-FORM LPS

TLR4 agonist

- LPS-B5 from *E. coli* O55:B5 **BEST**
- LPS-EB from *E. coli* O111:B4 **BEST**
- LPS-PG from *P. gingivalis* **BEST**

ROUGH (R)-FORM LPS

TLR4 agonist

- LPS-EK from *E. coli* K12
- LPS-SM from *Salmonella minnesota* R595

TLR4 antagonist

- LPS-RS from *Rhodobacter sphaeroides*

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TLR4 ANTAGONIST

LPS-RS from *Rhodobacter sphaeroides*

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To learn more agonistic and antagonistic effects of LPS on TLR4
Visit

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PRODUCTS

TLR4 ligands

PRODUCT	DESCRIPTION	QTY	CAT. CODE
TLR4 Agonist			
LPS-B5 Standard	Standard lipopolysaccharide from E. coli O55:B5	5 mg	tlrl-b5lps
LPS-B5 Ultrapure	Ultrapure lipopolysaccharide from E. coli O55:B5	5 mg	tlrl-pb5lps
LPS-EB Standard	Standard lipopolysaccharide from E. coli O111:B4	5 mg	tlrl-eb1ps
LPS-EB Ultrapure	Ultrapure lipopolysaccharide from E. coli O111:B4	5 x 10 ⁶ EU	tlrl-3pelps
LPS-EB Biotin	Biotinylated ultrapure LPS from E. coli O111:B4	500 µg	tlrl-lpsbiot
LPS-EB VacciGrade™	Preclinical grade ultrapure LPS from E. coli O111:B4	5 x 10 ⁶ EU	vac-3pelps
LPS-EK Standard	Standard lipopolysaccharide from E. coli K12	5 mg	tlrl-ek1ps
LPS-EK Ultrapure	Ultrapure lipopolysaccharide from E. coli K12	1 mg	tlrl-pek1ps
LPS-PG Standard	Standard lipopolysaccharide from P. gingivalis	1 mg	tlrl-pgl1ps
LPS-PG Ultrapure	Ultrapure lipopolysaccharide from P. gingivalis	1 mg	tlrl-ppg1ps
LPS-SM Ultrapure	Ultrapure lipopolysaccharide from S. minnesota	5 mg	tlrl-sm1ps
MPLA-SM	Monophosphoryl lipid A from S. minnesota	1 mg	tlrl-mpl1a
MPLA-SM VacciGrade™	Preclinical grade MPLA from S. minnesota	1 mg	vac-mpl1a
MPLAs	Synthetic monophosphoryl lipid A	1 mg	tlrl-mpl1s
MPLAs VacciGrade™	Preclinical grade synthetic MPLA	1 mg	vac-mpl1s
CRX-527	Synthetic lipid A analog	1 mg	tlrl-crx527
TLR4 Antagonist			
LPS-RS	Lipopolysaccharide from R. sphaeroides	5 mg	tlrl-rsl1ps
LPS-RS Ultrapure	Ultrapure lipopolysaccharide from R. sphaeroides	1 mg	tlrl-prsl1ps

Antibodies neutralizing TLR4-induced activation

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Anti-hTLR4-IgG	Human TLR4 neutralizing antibody - monoclonal mouse IgG1	100 µg	mabg-htlr4
PAb hTLR4	Human TLR4 neutralizing antibody - Polyclonal Rat IgG	200 µg	pab-hstlr4
Anti-hCD14-IgG	Neutralizing IgG monoclonal antibody to human CD14 (D3B8)	100 µg	mabg-hcd14

Synthetic inhibitors of TLR4-mediated signaling pathway

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Bay11-7082	I κ B- α inhibitor - NLRP3 inflammasome inhibitor	10 mg	tlrl-b82
BX795	TBK1/IKK ϵ inhibitor	5 mg	tlrl-bx7
CLI-095	TLR4 signaling inhibitor	1 mg	tlrl-cli95
MG-132	26S proteasome inhibitor - Autophagy activator	5 mg	tlrl-mg132
OxPAPC	TLR2 and TLR4 inhibitor	1 mg	tlrl-oxp1
Pepinh-MYD	MyD88 inhibitory peptide	2 mg	tlrl-pimyd
Pepinh-TRIF	TRIF inhibitory peptide	2 mg	tlrl-pitrif
Polymyxin B	LPS-induced TLR4 activation inhibitor	100 mg	tlrl-pmb
SB203580	p38 MAP Kinase Inhibitor	5 mg	tlrl-sb20
U0126	MEK1 and MEK2 Inhibitor	5 mg	tlrl-u0126

Related reporter cell lines

PRODUCT	DESCRIPTION	QTY	CAT. CODE
HEK-Blue™ hTLR4 Cells	Human TLR4 NF- κ B-SEAP reporter HEK293 cells	3-7 x 10 ⁶ cells	hkb-htlr4
HEK-Blue™ mTLR4 Cells	Murine TLR4- NF- κ B-SEAP reporter HEK293 cells	3-7 x 10 ⁶ cells	hkb-mtlr4
HEK-Dual™ mTLR4 Cells	Double readout murine TLR4 HEK293 cells	3-7 x 10 ⁶ cells	hkd-mtlr4ni
293/hTLR4a Cells	HEK 293 cells stably transfected with human TLR4a gene	3-7 x 10 ⁶ cells	293-htlr4a
293/mTLR4 Cells	HEK 293 cells stably transfected with murine TLR4 gene	3-7 x 10 ⁶ cells	293-mtlr4
293/hTLR4A-MD2-CD14 Cells	HEK 293 cells stably transfected with human TLR4a, MD2 and CD14 genes	3-7 x 10 ⁶ cells	293-htlr4md-2cd14
293/mTLR4-MD2-CD14 Cells	HEK 293 cells stably transfected with murine TLR4, MD2 and CD14 genes	3-7 x 10 ⁶ cells	293-mtlr4md-2cd14
293/hTLR4-HA Cells	HEK 293 cells stably transfected with HA-tagged human TLR4 gene	3-7 x 10 ⁶ cells	293-htlr4ha
THP1-Dual™ KO-TLR4 Cells	TLR4 knockout NF- κ B-SEAP and IRF-Lucia reporter human THP-1 monocytes	3-7 x 10 ⁶ cells	thpd-kotlr4
RAW-Dual™ KO-TLR4 Cells	TLR4 knockout IRF and MIP-2 (NF- κ B) reporter mouse macrophages	3-7 x 10 ⁶ cells	rawd-kotlr4
THP1-Dual™ KO-MyD Cells	MyD88 knockout NF- κ B-SEAP and IRF-Lucia reporter human THP-1 monocytes	3-7 x 10 ⁶ cells	thpd-komyd
RAW-Lucia™ ISG-KO-TRIF Cells	TRIF knockout IRF-Lucia reporter mouse macrophages	3-7 x 10 ⁶ cells	rawl-kotrif

Endotoxin detection

PRODUCT	DESCRIPTION	QTY	CAT. CODE
HEK-Blue™ LPS Detection Kit 2	Assay for the detection and quantification of biologically active LPS	1 kit	rep-lps2
HEK-Blue™ Selection	Antibiotics for maintenance of HEK-Blue™ Cells	10 x 1ml	hb-sel
HEK-Blue™ Endotoxin Standard	Standardized E. coli O55:B5 LPS	10 X 50 EU	rep-hbes-10
Normocin™	Antimicrobial reagent for preventing mycoplasma, bacteria and fungi	500 mg (10 x 1 ml) 1 g (1 x 20 ml)	ant-nr-1 ant-nr-2
HEK-Blue™ Detection	Cell culture medium for SEAP detection	5 pouches 10 pouches	hb-det2 hb-det3
QUANTI-Blue™ Solution	Alkaline phosphatase detection medium - Liquid form	5 ml (100X) 10 ml (100X) 20 ml (100X)	rep-qbs rep-qbs2 rep-qbs3

Selective Antibiotics

PRODUCT	DESCRIPTION	QTY	CAT. CODE
Blasticidin	Selective antibiotic for the bsr, bls, or BSD genes	100 mg (10 x 1 ml)	ant-bl-1
G418 (Geneticin)	Selective antibiotic for the neo gene	1 g (10 x 1 ml)	ant-gn-1
Hygromycin B Gold™	Selective antibiotic for the hph gene	1 g (10 x 1 ml)	ant-hg-1
Puromycin	Selective antibiotic for the pac gene	100 mg (10 x 1 ml)	ant-pr-1
Zeocin™	Selective antibiotic for the Sh ble gene	1 g (10 x 1 ml)	ant-zn-1

Visit
<https://www.invivogen.com/tlr4-ligands>
for products details

REFERENCES


1. Parusel, R. et al., 2017. An important question: Which LPS do you use?. *Virulence*, 8(8), 1890-1893.
2. Krukenberg, K. et al., 2015. Extracellular poly (ADP-ribose) is a pro-inflammatory signal for macrophages. *Chemistry & biology*, 22(4), 446-452.
3. Kiripolsky, J. et al., 2020. Activation of Myd88-Dependent TLRs Mediates Local and Systemic Inflammation in a Mouse Model of Primary Sjögren's Syndrome. *Frontiers in Immunology*, 10, 2963.
4. Evaristo-Mendonça, F. et al., 2019. Preconditioning of Rat Bone Marrow-Derived Mesenchymal Stromal Cells with Toll-Like Receptor Agonists. *Stem cells international*, 2019.
5. Bi, D. et al., 2018. Recognition of lipopolysaccharide and activation of NF- κ B by cytosolic sensor NOD1 in teleost fish. *Frontiers in immunology*, 9, 1413.
6. Cao, Z. et al., 2018. Tissue expression of Toll-like receptors 2, 3, 4 and 7 in swine in response to the Shimen strain of classical swine fever virus. *Molecular Medicine Reports*, 17(5), 7122-7130.
7. Vogel, D. Y. et al., 2014. Human macrophage polarization in vitro: maturation and activation methods compared. *Immunobiology*, 219(9), 695-703.
8. Nagala, M. et al., 2018. Expression of Siglec-E alters the proteome of lipopolysaccharide (LPS)-activated macrophages but does not affect LPS-driven cytokine production or toll-like receptor 4 endocytosis. *Frontiers in immunology*, 8, 1926.
9. Seppola, M. et al., 2015. Ultrapure LPS induces inflammatory and antibacterial responses attenuated in vitro by exogenous sera in Atlantic cod and Atlantic salmon. *Fish & shellfish immunology*, 44(1), 66-78.
10. Rialdi, A. et al., 2016. Topoisomerase 1 inhibition suppresses inflammatory genes and protects from death by inflammation. *Science*, 352(6289).
11. Guida, C. et al., 2015. A novel inflammatory pathway mediating rapid hepcidin-independent hypoferrremia. *Blood, The Journal of the American Society of Hematology*, 125(14), 2265-2275.
12. Herbst-Kralovetz, M. M. et al., 2013. Lack of norovirus replication and histo-blood group antigen expression in 3-dimensional intestinal epithelial cells. *Emerging infectious diseases*, 19(3), 431.
13. Park, J. et al., 2017. Fucoidan inhibits LPS-induced inflammation in vitro and during the acute response in vivo. *International immunopharmacology*, 43, 91-98.
14. Chang, P. C. et al., 2013. Irradiation by light-emitting diode light as an adjunct to facilitate healing of experimental periodontitis in vivo. *Journal of periodontal research*, 48(2), 135-143.
15. Santegoets, K. C. et al., 2016. Impaired Porphyromonas gingivalis-induced tumor necrosis factor production by dendritic cells typifies patients with rheumatoid arthritis. *Arthritis & Rheumatology*, 68(4), 795-804.
16. Gözl, L. et al., 2015. Hypoxia and P. gingivalis synergistically induce HIF-1 and NF- κ B activation in PDL cells and periodontal diseases. *Mediators of inflammation*, 2015.
17. Feng, Z. et al., 2018. METTL 3 regulates alternative splicing of MyD88 upon the lipopolysaccharide-induced inflammatory response in human dental pulp cells. *Journal of cellular and molecular medicine*, 22(5), 2558-2568.
18. Chen, Z. et al., 2015. IL-1R/TLR2 through MyD88 divergently modulates osteoclastogenesis through regulation of nuclear factor of activated T cells c1 (NFATc1) and B lymphocyte-induced maturation protein-1 (Blimp1). *Journal of Biological Chemistry*, 290(50), 30163-30174.
19. Saz-Leal, P. et al., 2018. Targeting SHIP-1 in myeloid cells enhances trained immunity and boosts response to infection. *Cell reports*, 25(5), 1118-1126.
20. Kang, S. S. et al., 2018. Staphylococcal LTA antagonizes the B cell-mitogenic potential of LPS. *Scientific reports*, 8(1), 1-12.

REFERENCES (CONT'D)

21. McAlees, J. W. et al., 2015. Distinct Tlr4-expressing cell compartments control neutrophilic and eosinophilic airway inflammation. *Mucosal immunology*, 8(4), 863-873.
22. Rosa, A. et al., 2015. HIV-1 Nef promotes infection by excluding SERINC5 from virion incorporation. *Nature*, 526(7572), 212-217.
23. Carvelli, J. et al., 2020. Association of COVID-19 inflammation with activation of the C5a-C5aR1 axis. *Nature*, 1-9.
24. Kroes, M. M. et al., 2019. Activation of human NK cells by *Bordetella pertussis* requires inflammasome activation in macrophages. *Frontiers in Immunology*, 10, 2030.
25. Ferdous, F. et al., 2008. Broiler chick thrombocyte response to lipopolysaccharide. *Poultry science*, 87(1), 61-63.
26. Berthet, J. et al., 2012. Human platelets can discriminate between various bacterial LPS isoforms via TLR4 signaling and differential cytokine secretion. *Clinical Immunology*, 145(3), 189-200.
27. Brykczynska, U. et al., 2020. Distinct Transcriptional Responses across Tissue-Resident Macrophages to Short-Term and Long-Term Metabolic Challenge. *Cell Reports*, 30(5), 1627-1643.
28. Ravindran, R. et al., 2016. The amino acid sensor GCN2 controls gut inflammation by inhibiting inflammasome activation. *Nature*, 531(7595), 523-527.
29. Krikun, G. et al., 2012. Lipopolysaccharide appears to activate human endometrial endothelial cells through TLR-4-dependent and TLR-4-independent mechanisms. *American journal of reproductive immunology*, 68(3), 233-237.
30. Schürmann, M. et al., 2020. Stem Cell-Induced Inflammation in Cholesteatoma Is Inhibited by the TLR4 Antagonist LPS-RS. *Cells*, 9(1), 199.
31. Abdelmageed, M. E. et al., 2016. LPS-RS attenuation of lipopolysaccharide-induced acute lung injury involves NF- κ B inhibition. *Canadian Journal of Physiology and Pharmacology*, 94(2), 140-146.
32. Kawakita, F. et al., 2017. Effects of Toll-like receptor 4 antagonists against cerebral vasospasm after experimental subarachnoid hemorrhage in mice. *Molecular neurobiology*, 54(8), 6624-6633.
33. Tan, Y. et al., 2015. Mechanisms of Toll-like receptor 4 endocytosis reveal a common immune-evasion strategy used by pathogenic and commensal bacteria. *Immunity*, 43(5), 909-922.
34. Nguyen, Q. T. et al., 2020. E. coli-Produced Monophosphoryl Lipid a Significantly Enhances Protective Immunity of Pandemic H1N1 Vaccine. *Vaccines*, 8(2), 306.
35. Han, X. et al., 2019. In situ thermal ablation of tumors in combination with nano-adjuvant and immune checkpoint blockade to inhibit cancer metastasis and recurrence. *Biomaterials*, 224, 119490.
36. Cheng, Q. et al., 2019. Sequential conditioning-stimulation reveals distinct gene-and stimulus-specific effects of Type I and II IFN on human macrophage functions. *Scientific reports*, 9(1), 1-14.
37. Hellmann, I. et al., 2018. Novel antibody drug conjugates targeting tumor-associated receptor tyrosine kinase ROR2 by functional screening of fully human antibody libraries using Transpo-mAb display on progenitor B cells. *Frontiers in Immunology*, 9, 2490.
38. Albershardt, T. C. et al., 2020. Intratumoral immune activation with TLR4 agonist synergizes with effector T cells to eradicate established murine tumors. *npj Vaccines*, 5(1), 1-9.
39. Yang, J. et al., 2020. Melatonin pretreatment alleviates renal ischemia-reperfusion injury by promoting autophagic flux via TLR4/MyD88/MEK/ERK/mTORC1 signaling. *The FASEB Journal*, 34(9), 12324-12337.
40. Li, X. D. et al., 2020. Non-canonical Activation of Human Group 2 Innate Lymphoid Cells by TLR4 Signaling. *bioRxiv* [pre-print].

LPS

INVIVOGEN'S PRODUCTS FOCUS

www.invivogen.com

InvivoGen Asia

Unit 106; 1F, 15W Phase 3
Hong Kong Science Park, Pak Shek Kok,
Hong Kong, China

Toll-free: +852-36223480

Tel: +852-36223480

Fax: +852-36223483

Email: info.hk@invivogen.com