

## ■金属スカベンジャー対応表

Scavenger	Ag	Ca	Cd	Co	Cr	Cs	Cu	Fe	Hg	Ir	La	Li	Mg	Ni	Os	Pb	Pd (II)	Pd (0)	Pt	Rh (I)	Rh (II)	Rh (III)	Ru (II)	Sc	Sn	W	Zn
Si-Thiol	■						□		■	□					■	□	■	■		□	□	□	□	□	□		
Si-Thiourea	□						□	□							□		■	■		□	□	□		□	□		
Si-Cysteine		□	■		□	□	□	■		■	□		□		■		□	□	□	□	□		■	■	■		□
Si-DMT			□	□			□	□		■				■	■		■	■	■	■	■	■	■	□			□
Si-Amine			■	□	■		□	□	□							□	□		■	■	■					□	□
Si-Diamine			□	□	■		□	□	□					□		□	■	■	■					□		□	□
Si-Triamine				□	■		□	□						□		■	■	■	■					□		□	□
Si-Imidazole			■	■	□		■	■						■	■		□	□		□	□					■	■
Si-TAAcOH				■	□	□		□						■	■		□	□		□	□			■	□		
Si-TAAcONa		■	■		□	■	■	■		■	■	■	■	■	■		□	□		□	□	■	□	■	■		□
Si-TBD		-		■	■	-	■			-	-	-	-		-		□							-		-	

■ Preferred scavengers □ Scavenges - 該当データなし