



NACALAI TESQUE, INC.
Nijo Karasuma, Nakagyo-ku, Kyoto 604-0855, JAPAN

Product No : 06881 Grade :
Product Name : Magnesium Sulfate, Anhydrous

MF : MgSO₄
FW : 120.37

Test item	Unit	Specified value	Test method
Purity (after ignition)	%	min. 99.0	
Solubility in water		Clear to faint turbidity	
Loss on ignition	%	max. 4	
Chloride (Cl)	%	max. 0.02	
Lead (Pb)	%	max. 0.001	
Arsenic (As)	ppm	max. 4	